



**Medical Footcare**  
Foot Specialists, Since 1885

**New!**  
**1stRaythotics™  
Preforms®**  
New Treatment  
Options for Turf  
Toe, Hallux Limitus,  
Rigidus, Arthritis,  
Post-Op, More  
#4870 **PAGE 6**

**New!**  
**1/2" Lateral Sole  
Wedge Insoles**  
Limit Supination,  
Relieve Peroneal  
Tendon, Knee  
Pain, More  
#2301 **PAGE 11**

**New!**  
**PressureOFF® Customizable Insoles**  
Offload Plantar Pathologies  
in Everyday Footwear  
#2408 **PAGE 4**

**New!**  
**Visco-GEL® Customizables®**  
The First & Only Ortho-digital Devices  
You Personalize For Each Patient  
#1030-CTZ **PAGE 28**

**New!**  
**Just For Toenails®**  
Enriched Polish Improves the  
Appearance of Problem Nails  
#M905 **PAGE 60**

**New!**  
**Springthotics™ Spring Plates**  
**PAGE 10**

**New!**  
**Halluxthotics™  
Turf Toe Arch Supports**  
**PAGE 9**

**New!**  
**FungaSoap® CBD Relief  
Formula** **PAGE 62**

**Unique Foot Treatments to Help You,  
Your Patients & Your Practice**  
Dispense • Prescribe • Patient Direct Purchase

**30 New Products...  
See Inside!**

Volume 22M To order or get more information,  
Call: 1-800-424-5561 • Visit: [medical.pedifix.com](http://medical.pedifix.com)  
Email: [info@pedifix.com](mailto:info@pedifix.com) • Fax: 1-800-431-7801  
Distributors Nationwide

## Offloading & Wound Care Innovation



**New!**

## PressureOFF® Customizable Offloading Insoles

*Help More Patients  
Relieve Plantar Pathologies  
with the  
First 'Removable Pegs'  
Offloading Insoles  
for Everyday Footwear*

- Proven to Reduce Pressure, Friction & Pain
- Likely to Save Limbs & Lives
- In the Shoes Your Patients Will Wear
- For Higher Offloading Compliance & Better Outcomes
- Exclusive, Patented Designs



*"I'll use these to treat a variety of conditions more effectively, and with higher patient acceptance than other offloading options"*

*John Powers, DPM*

*"Many of my patients benefit from offloading, and these insoles make it easy, economical and convenient to relieve plantar pathologies"*

*Tim Short, DPM*

## The Removable Pegs Insole Design You Know, Now for Regular Shoes

PressureOFF® Insoles enable targeted offloading of plantar pathologies in your patient's footwear.

The proven 'removable pegs' design — previously available only for post-op and wound-care shoes — helps prevent, relieve and promote healing of dozens of common conditions.

### New Treatment Options for:

Calluses

Corns

Blisters

Bruises

Warts

Prominences

Bone Spurs

Lesions

Wounds

IPKs

Porokeratosis

Metatarsalgia

Sesamoiditis

Pressure Points

Painful Spots

Fat Pad Atrophy

Sensitive Areas

Friction Zones

Pre-Ulcer 'Hot Spots'

Diabetic Foot Ulcers

Ulcers in Remission

Surgical Sites

More

- Removable Plantar Pegs Unload Targeted Areas
- Offload Sore Spots to Relieve Pain, Reduce Friction
- Minimize Skin Breakdown Due to Repetitive Micro Trauma
- Distribute Pressure Away from Sensitive Areas
- Prevent, Relieve, Heal

If you're treating or preventing plantar pathologies and diabetic ulcers, get your patients on PressureOFF® Customizable Offloading Insoles. The easy, convenient, discrete offloading modality they'll love — and hardly notice.



See Your Choices →

More information: [www.pedifix.com/t-POIOffloading.aspx](http://www.pedifix.com/t-POIOffloading.aspx)

## A New Offloading Modality Patients Comply With

### Introductory PressureOFF® Offloading Systems

Choose Insole Characteristics That Match Your Patient's Needs

**Things to consider:**

- Insole Thickness
- Offloading Heights
- Material Layers
- Shoe Fit



**Sizing Guide — Size-up for wider foot coverage**

Specify Right or Left

- Extra Large** fits shoe sizes 13-15 (+/-31.5cm),
- Large** fits 10-12 (+/-28.6cm),
- Medium** fits 6-9 (+/- 25.7cm).

**Important** – Scissor-Trim Precisely to Match Foot, Footbed or Original Shoe Insole Size & Shape



**PressureOFF®**

**High Protection**  
Offload Factor 6/9

9.0mm Total Cushioning

- 3.0mm Poron®
- 3.0mm EVA +
- 3.0mm EVA Pegs


6.0mm Offloading

**Design Distinctions\***

- Most durable top layer
- Tallest removable peg height
- Maximum offload
- The removable Peg Pattern excludes an area just distal the heel.

Specify Size, Right or Left

#2402 1/pk  
POI System 9



**PressureOFF®**

**High Protection**  
Offload Factor 6/9

**Wide Shape**

9.0mm Total Cushioning

- 3.0mm Poron®
- 3.0mm EVA +
- 3.0mm EVA Pegs


6.0mm Offloading

**Design Distinctions\***

- Most durable top layer
- Tallest removable peg height, maximum offload
- WIDE SHAPE for wide & deformed feet
- The removable Peg Pattern excludes an area just distal the heel.

Specify Size, Right or Left

#2404 1/pk  
POI System 9 Wide



### Supplemental Top Covers

Extra Foot Protection  
Reinforce or Refresh Top Insole Layer



Quick, easy, self-sticking application



**2.0mm Plastazote® Supplemental Top Cover**

Self-Adhering, One Size, Trim to Fit

Right (#2415-XLR)  
Left (#2415-XLL)



**2.0mm Poron® Supplemental Top Cover**

Self-Adhering, One Size, Trim to Fit

Right (#2416-XLR)  
Left (#2416-XLL)

### PegSaver™ Base



**1.5mm Stabilizer Base**

Included with every PressureOFF® insole, these self-adhering bottom layers are optional to apply, but highly recommended to extend insole life, slow peg compression and migration.

One, Trim-to-Fit Size  
Standard Shape, #2410, Left or Right  
Wide Shape, #2411, Left or Right



**PressureOFF®**

**Offload Factor 3/7**

**Retachable™ Pegs**

Poron®/Plastazote®

**7.0mm Total Cushioning**

1.5mm Poron®

1.5mm Plastazote®

1.0mm MicroLoop™

3.0mm Poron® Pegs

**3.0mm Offloading**

**Design Distinctions\***

- MicroLoop™ Pegs — Adjust, Replace Easily as Pathologies Change
- Dual Top Layer Options — Choose Poron® or Plastazote®
- Removable Peg Pattern extends nearly to insole perimeter.
- Poron® Pegs may prove most durable.

Specify Size, Right or Left

**#2406 1/pk**

POI System 7



**PressureOFF®**

**Offload Factor 3/7**

**Retachable™ Pegs / Soft Top**

Plastazote®/Poron®

**7.0mm Total Cushioning**

1.5mm Plastazote®

1.5mm Poron®

1.0mm MicroLoop™

3.0mm Poron® Pegs

**3.0mm Offloading**

**Design Distinctions\***

- MicroLoop™ Pegs Adjust, Replace Easily as Pathologies Change
- Dual Top Layer Options — Choose Plastazote® or Poron®
- Removable Peg Pattern extends nearly to insole perimeter.
- Poron® Pegs may prove most durable.

Specify Size, Right or Left

**#2408 1/pk**

POI System 7 Soft

**PressureOFF®**

**Soft Top**

**Offload Factor 4/6**

**6.0mm Total Cushioning**

2.0mm Plastazote®

2.0mm EAV +

2.0mm Poron® Pegs

**4.0mm Offloading**

**Design Distinctions\***


- Thickest Plastazote® Top
- Thinnest overall design
- Removable Peg Pattern does not extend to the insole perimeter. Poron® Pegs may prove most durable. Super soft — watch for ring edema, edge abrasion.

One, Extra Large, Trimmable Size Fits Up to Men's 15, 31.2cm long.

Specify Right or Left

**#2400 1/pk**

POI System 6 Soft



**Dual Offloading Advantages!**

**1.** Customize PressureOFF® bottoms by removing pegs under specific areas.

**2.** Customize PressureOFF® tops with self-adhering Metatarsal, Arch, 'Dancer's' and 'U'-shaped Pads. See Gel, Felt and Foam options — pages 15, 17, 37–39.



\*Please Note — these new designs are being refined in response to clinician and patient feedback, so please share yours. During their introduction and evolution — amidst worldwide shortages of Foam supplies as this catalog goes to press — materials and thicknesses may vary from those specified here. Thank you for your understanding as we improve these products with your input. Help us make these better... email info@pedifix.com

More information: [www.pedifix.com/t-POIOffloading.aspx](http://www.pedifix.com/t-POIOffloading.aspx)

New Hallux Treatment Options



**New!**

**1stRaythotics™  
Preforms®**

**Integrated Morton's Extensions  
& Mid-Foot Support**

**Stabilize, Limit First Ray Motion  
to Relieve Big Toe Joint Pain**

**Carbon Fiber Thin for Easy Shoe Fit**

Combining Morton's Extensions with Integrated Medial Support and High-Performance Carbon Fiber, these unique orthotics provide comfortable relief and promote healing of various hallux and forefoot ailments.

1stRaythotics™ offer midfoot support while they stabilize and limit first metatarsal joint range of motion to relieve pain. Thinner and lighter than steel, our exclusive, German-made Carbon Fiber material offers just the right amount of hallux flexibility and arch support, and fit in most shoes.

- **Limit Joint Motion**
- **Reduce Shoe Flexibility**
- **Enhance Foot Function**
- **Speed Healing**

**Indications**

- Turf Toe
- Hallux Limitus
- Hallux Rigidus
- Hallux Valgus
- Hallux Arthritis
- Sesamoiditis
- Plantar Plate
- Metatarsalgia
- Morton's Neuroma
- Hammer Toes
- Flatfoot
- Diabetic Hot Spots
- Post-Op Support
- More



**Never dispense  
another flat  
TT plate again!**



Thin, Comfortable, Flexible Support

More information: [www.pedifix.com/t-1stRaythotics.aspx](http://www.pedifix.com/t-1stRaythotics.aspx)



**Morton's Extension  
Bilateral**  
#4870



**Morton's Extension Left Foot,  
Right Foot Neutral**  
#4871



**Morton's Extension Right Foot,  
Left Foot Neutral**  
#4872

**Neuroma Relief!**



**Reverse Morton's Extension  
Bilateral**  
#4875



**Metatarsal Offloader  
Bilateral**  
#4874



**Neutral Bilateral**  
#4873

**1pair/pk** 10 Sizes — M4/W5, M5/W6, M6/W7, M7/W8, M8/W9, M9/W10, M10/W11, M10 ½, M11, M12, M13, M14/M15



**New!**

**1stRayometer™ First Ray Measuring Device**

*Assess & quantify first ray mobility  
and range of motion consistently*

- Anatomically curved, easy, fast, precise, consistent measurement of dorsi and plantar flexion in millimeters.
- Measure, demonstrate, document dysfunction
- Substantiate treatment decisions
- Helpful for Hallux Limitus/Rigidus, Turf Toe, Arthritis, Sesamoiditis and other conditions.



**1/pk Item#**  
#9000

**Stabilize, Limit First Toe Joint Motion to Relieve Pain**



**New!**

**1stRaythotics™  
Hallux Trainer  
Sports Insoles**

**'Morton's Extension' Support Inside Orthotics Relieves Big Toe Pain**

**Treat Turf Toe, Hallux Limitus, Hallux Rigidus, Arthritis, Sesamoiditis, Post-Op, More**

For many patients with first MTPJ conditions, this is the most comfortable solution. The 'hidden' Morton's Extension inside these orthotics stabilize and limit first toe joint motion to ease pain. Sturdy medial arch support, cushioned top overs, soft metatarsal pads and deep heel cups redistribute weight, pressure and biomechanical stresses to relieve pain and promote healing of hallux conditions. These replace existing insoles in footwear. See sizing information below.

**1 pr/pk Item#**  
**#4876**



U.S. Shoe Size/ Euro Size	Insole Length	U.S. Shoe Size/ Euro Size	Insole Length
W4/Euro 35	9"	W11/M10/Euro 42	11"
W5/Euro 36	9 1/4"	W12/M11/Euro 43	11"
W6/M5/Euro 37	9 3/4"	M12/Euro 44	11 1/2"
W7/M6/Euro 38	10"	M13/Euro 45	11 3/4"
W8/M7/Euro 39	10 1/4"	M14/Euro 46	12"
W9/M8/Euro 40	10 1/2"	M15/Euro 47	12 1/4"
W10/M9/Euro 41	10 3/4"	M16/Euro 48	12 1/2"



**New!**

**Halluxthotics™  
GXT Fiberglass  
Turf Toe  
Arch Supports**

*'Morton's Extension' Stabilizes,  
Limits Big Toe Joint Motion to  
Relieve Pain Much More  
Comfortably Than Flat TT Plates.*



Stiffen shoes to reduce their flexibility and big toe joint movement to ease Turf Toe, Hallux Limitus, Rigidus, Arthritis, and other 1st toe specific conditions. Molded arch for support and comfort. Very thin, fit easily in most footwear. Sold in pairs. See sizing information below.

**1pr/pk Item#**  
**#4868**

U.S. Shoe Size/ Euro Size	Insole Length	U.S. Shoe Size/ Euro Size	Insole Length
W5/M4/Euro 35/36	8 1/2"	W10-11/M9-10/Euro 41/42	9 5/8"
W6-7/M5-6/Euro 37/38	8 7/8"	M11-12//Euro 43/44	10 1/4"
W8-9/M7-8/Euro 39/40	9 1/4"	M13/M14/Euro 45/46	10 3/4"



**New!**

**Vasyli Dananberg Insoles**

*Improve Big Toe Movement to Soothe Functional Hallux Limitus  
First with removable wedges & 3-plane motion control*

Unique orthotics with adjustable 'tabs' under the great toe joint to personalize range of motion, improve first ray function and relieve pain. Customize these to comfort and get them back to their favorite activities faster. Solid arch support, heel cupping shock absorption and personalized control of the big toe joint also help relieve Bunion Pain, Plantar Fasciitis, Posterior Tibialis Tendonitis, even low back pain.

**1pr/pk Item#**  
**#4850**

**Sizing Guide**

XSmall fits Women's 5-6  
Small Fits Women's 7-8, Men's 6-7  
Medium fits Women's 9-10, Men's 8-9  
Large fits Women's 11-12, Men's 10-11  
XLarge fits Men's 12-13

**New!**

**Halluxthotics™ XTC Carbon Fiber  
Graphite Turf Toe Relief Insole**

*Flat Insert Limits Big Toe Joint  
Motion That Causes Pain*

'Morton's Extension' supports, stabilizes, reduces big toe joint range of motion to ease pain. Stiffens shoes, reduces flexibility under the big toe specifically. Fit easily in any footwear style. Sold singly, not pairs. Flip for Right or Left foot use. Just two sizes — be sure to choose the length that is shorter than your shoe insole. Position these under an insole or orthotic. Small Size is 8". Large Size is 10 1/4".



**1/pk Item#**  
**#4862**

**New!**

**SOFT Shoe Insoles  
Cushioned Comfort and  
as Top Covers**

Position these in thin-soled shoes for extra padding, or on top of the shoe stiffening insoles shown here. Scissor trimmable. Specify size.



**1pr/pk Item#**  
**#2330**

U.S. Shoe Size/ Euro Size	U.S. Shoe Size/ Euro Size	U.S. Shoe Size/ Euro Size
W6/Euro 36	W11/M8/Euro 41	M12/Euro 45
W7/Euro 37	W12/M9/Euro 42	M13/Euro 46
W8/Euro 38	M10/Euro 43	M14/Euro 47
W9/Euro 39	M11/Euro 44	M15/Euro 48
W10/M7/Euro 40		

**New!**

**GXT Fiberglass Full  
Foot Insoles — Flat**

*Maximum Rigidity to Limit Shoe Flexibility, Reduce Stress on Toes,  
Forefoot & Midfoot Conditions*

Rigid Inserts make shoes really stiff — almost inflexible — limiting joint motion. Very thin, to fit easily in most footwear. No arch support or toe spring contour. Not recommended for athletic use — too rigid. Medium width only. See sizing information below. Choose an Insole size smaller than the shoe size.



**1pr/pk Item#**  
**#4866**

U.S. Shoe Size/ Euro Size	Insole Length	U.S. Shoe Size/ Euro Size	Insole Length
W6/M5/Euro 37	9 3/4"	W11/M10/Euro 42	11"
W7/M6/Euro 38	10"	W12/M11/Euro 43	11 1/4"
W8/M7/Euro 39	10 3/8"	M12/Euro 44	11 1/2"
W9/M8/Euro 40	10 1/2"	M13/Euro 45	12"
W10/M9/Euro 41	10 3/4"		



Relieve, prevent, promote healing of these common conditions



- Big Toe Joint Pain
- Turf Toe Injury
- Hallux Limitus
- Hallux Rigidus
- Toe & Forefoot Arthritis
- Forefoot Metatarsal Pain
- Forefoot Calluses
- Forefoot Capsulitis, Bursitis
- Foot Sprains, Fractures
- Metatarsalgia
- Midfoot Conditions
- Morton's Neuroma
- Freiberg's Infraction
- Arch & Heel Pain
- Sesamoiditis/Injury
- Cuboid Syndrome
- Achilles Tendinitis
- Plantar Plate Injuries
- Diabetic Foot Ulcers
- Trans Metatarsal Amputations
- Gait Abnormalities
- Plantar Puncture Protection
- Post-Surgical Stabilization
- More



Relieve pain from hyperextension injuries

**New!**

## Springthotics™ Carbon Fiber Graphite Spring Plates

Stiffen Shoes to Reduce Pain, Propel You Forward — Comfortably

Relieve, support, promote healing of big toe, forefoot and midfoot conditions, sprains or fractures. Exclusive material and design reduce shoe flexibility to stabilize joint motion that causes pain. Rocker Sole and Toe Spring shape enhance ambulation and energy transfer through the gait cycle to make walking easier. So thin, they'll fit in most shoes, and under orthotics or insoles.

**Thinner,  
Lighter,  
Stronger  
Than Steel!**



### XLT Carbon Fiber Spring Plate Graphite Insole (Semi-Rigid)

Sold singly, for Right or Left Foot — be sure to specify. If you need a pair, order 2 (R & L). See IMPORTANT sizing information below\*.

1/pk	Item#		
	<b>#4864</b>		
Size	Insert Length	Size	Insert Length
W5	9 1/8"	M7	10"
W6	9 3/8"	M8	10 1/4"
W7	9 11/16"	M9	10 3/4"
W8	10"	M10	10 7/8"
W9	10 1/4"	M11	11 1/2"
W10	10 9/16"	M12	11 3/4"
		M13	12 1/16"



'Rocker Sole' with Toe Spring Propels You Forward



**"Very effective in preventing toe pain. I can now walk and exercise again!"**

★★★★★

**\*Important sizing information:** These orthotics are available in Medium width only. They will not fit into narrow footwear, and they cannot be trimmed. Instead of choosing by shoe sizes, which vary, measure the length of the insole in your favorite shoes — from the heel to the tip of the toe. Use that to pick your full-length insole size, making certain it is shorter than your existing insole so it will fit in your shoes. If you wear a half-size, or are unsure, size down to the next smaller size Insert — they need not fill your footwear to be effective. Tips: It is important not to wear these in shoes that have very flexible soles, as the lack of rigidity will stress, weaken or crack some Inserts. \*For maximum comfort, these are intended for placement under an insole, shoe liner or orthotics. Most types of insole will be effective as a top cover.



### XAT Carbon Fiber Spring Plate Graphite Insole (Rigid)

Sold in Pairs. Slightly stiffer than XLT, to further reduce shoe flexibility, decrease range of motion that hurts. See IMPORTANT sizing information below\*.

**1 pr/pk Item#**  
**#4863**

U.S. Shoe Size/ Euro Size	Insole Length	U.S. Shoe Size/ Euro Size	Insole Length
W5/M4/Euro 36	9 1/2"	M11/Euro 43	11 1/4"
W6/M5/Euro 37	9 3/4"	M12/Euro 44	11 1/2"
W7/M6/Euro 38	10"	M13/Euro 45	12"
W8/M7/Euro 39	10"		
W9/M8/Euro 40	10 1/2"		
W10/M9/Euro 41	10 3/4"		
W11/M10/Euro 42	11"		



Coils Energy When Your Foot Lands, Springs as You Toe-Off





**New!**

**Lateral Sole Wedge Insoles**

*Limit Excess Supination Symptoms, Peroneal Tendon, Foot & Knee Pains*

These beveled, 'sloped' insoles create a valgus wedge under the entire lateral aspect of the foot. They help control supination, ease Peroneal Tendon strain, help stabilize ankles, decrease uneven shoe wear, and reduce medial knee pain.

**1/2" thick** varus/valgus correction. Real Leather top, leatherette bottom, durable Polyurethane core.

Choose Sizes:

Small: Euro 35-36, fits Women's Shoe Sizes 5, Men's 4 (9 3/8" long)

Medium: Euro 37-38, fits Women's Shoe Sizes 6-7, Men's 5-6 (9 3/4" long)

Large: Euro 39-40, fits Women's Shoe Sizes 8-9, Men's 7-8 (10 1/8" long)

XLarge: Euro 41-42, fits Women's Shoe Sizes 10-11, Men's 9-10 (10 1/2" long)

XXLarge: Euro 43-44, fits Women's 12-13, Men's 11-12 (11" long)

**1pr/pk Item#**  
**#2301**

**1/4" thick**, durable Nylon cover, Poron® middle layer, sturdy Polyurethane base.


Choose Sizes:

Small/Medium: Fits Women's Shoe Sizes 5-9, Men's 6-9 (11" long x 3 1/2" wide)

Large/XL: Fits Women's Shoe Sizes 9+, Men's 10-12 (12" long x 3 3/4" wide)

XXLarge: Fits Men's 13+ (12.5" long x 4" wide)

**1pr/pk Item#**  
**#2300**



# SoreSpot<sup>®</sup> Silicone Tape



**Waterproof!**  
Stays on after  
shower or bath

- Reusable
- Hypoallergenic
- Latex-free

## Skin Saver!


Foot Care, Skin Care,  
Wound Care,  
Scar Care  
— Anywhere!

**Unique, Self-Adhering  
Protection for  
Sensitive Skin**

- Shields Blisters, Corns & Irritations
- Reduces Friction & Skin Damage
- Gently Secures Bandages & Dressings
- Helps Improve the Appearance of Scars
- More!

**Safe & Easy Release for Painless  
Removal Without Skin Damage**

### Scars



Before

After



**“My patients call it ‘Miracle Tape’”**

*Dr. Michele Kurlanski, Lighthouse Foot & Ankle*

**SoreSpot® Silicone Tape is a ground-breaking advance in skin protection, with multiple benefits and uses**

First, this unique skin protectant minimizes shear and friction to relieve and prevent blisters, corns, calluses, irritations and 'hot spots', for instant pain relief anywhere on the body. It is helpful for young and old, sedentary and active customers, including runners and other athletes.

Second, silicone is the most clinically-proven solution for improving the appearance of scars, old or new. SoreSpot® Silicone Tape makes it easy and economical to benefit from its noticeably visible effects.

Third, silicone creates an occlusive barrier over skin that enhances the penetration of wart topicals for faster treatment results.

Fourth and perhaps most valuable, unlike paper and cloth tapes, silicone releases so gently that it prevents tears and damage to fragile skin upon removal, making it ideal for positioning bandages on elderly and at-risk patients, including those with PVD, diabetes, and other concerns.

Soft, flexible and conforming, SoreSpot® Silicone Tape wraps easily around toes, protects 'pump bumps' and bunions, too. The tape is waterproof, so it stays in place through showers and baths, even swimming. Hypo-allergenic and latex-free, it won't irritate skin or cause reactions.

So soothing, users report success also on burns, poison ivy, bug bites and other skin conditions, too.

SoreSpot® Silicone Tape is now available in 1- and 2-inch widths, 1.5- and 5-yard lengths, in single rolls or 12-roll 'clinic boxes'. Discover its many uses and benefits for yourself and your patients.

Link to SoreSpot® Silicone Tape Video on YouTube here:  
<https://www.youtube.com/watch?v=yq3sF34qZX8>

1 Roll	Item#
1 in. x 1.5 yd. OTC Pack	#P3095-115
1 in. x 1.5 yd. Medical Pack	#3095-115
1 in. x 5 yd. Medical Pack	#3095-150
2 in. x 1.5 yd. Medical Pack	#3095-215
2 in. x 5 yd. Medical Pack	#3095-250

12 Roll Clinic Pack	Item#
1 in. x 1.5 yd.	#3095-115-12*
1 in. x 5 yd.	#3095-150-12*
2 in. x 1.5 yd.	#3095-215-12*
2 in. x 5 yd.	#3095-250-12*

\*Not packaged for individual resale

**Medical Polybag Pack**



**OTC Blister Pack**  
 (1 in. x 1.5 yd.)

**12-Packs**





## Exclusive Skin Protection Options



- Sturdy
- Tear-Resistant
- Self-Adhering
- Stays in Place



Help Improve the Appearance of Scars!


*“Silicone has the best evidence of all treatments for scars”*

Mustoe et al, 2002

### Visco-GEL® Enhanced Silicon™ Sheets & SkinStrips™

*Protects Skin, Reduces Scars*

Two benefits and applications! Our special silicon — enriched with SmartGEL® additives — protects skin from shear and minimizes friction to shield calluses and other sore spots. Also helps smooth, soften and improve the texture, color and thickness of scars. Self-adhering. Washable and reusable for repeated applications. Soft, pliable, comfortable. Fabric top cover increases durability, and further reduces friction. Cut to any shape!




	Item#
2 1/2" x 3" x 2.0mm Sheet, Beige Top Cover	10/pk #805
2 1/2" x 3" x 2.0mm Sheet, Clear Top Cover	10/pk #806
2 1/2" x 3" x 2.0mm Sheet, Beige Top Cover	2/pk #3075
2 1/2" x 3" x 2.0mm Sheet, Clear Top Cover	2/pk #3076
4" x 6" x 2.0mm Sheet, Beige Top Cover	1/pk #3080
4" x 6" x 2.0mm Sheet, Clear Top Cover	1/pk #3081
1" x 8" x 1.0mm SkinStrip, Beige Top Cover	1/pk #3091
2" x 8" x 1.0mm SkinStrip, Beige Top Cover	1/pk #3092




#3091

*SkinStrips™ are Ultra Thin, just 1.0mm!*



#3092



#3075 (#805)



#3076 (#806)



#3080



#3081

**Relieve. Remove. Wash. Reuse.**

**Pedi-GEL® Reuseable Self-Sticking Shoe & Foot Pads**

*Cushion and Protect to Relieve Pain*

This unique new material adheres directly to skin or to the insides of shoes— whichever is preferred—to provide extra foot cushioning and support. Self-sticking Gel is washable and reusable, providing long-lasting relief from foot pain or shoe pressure. Choose a pre-made shape to relieve a specific discomfort, or make your own from flat sheets.

**Ball-of-Foot Pads**

- Item #8201, 2/package

**Arch Pads**

- Item #8202, 2/package

**Callus Pads**

- Item #8203, 2/package

**Dancer's Pads**

- Item #8204, 2/package

**Comfort Gel Sheets**

- Item #8205, 2/package

**Corn Cushions**

- Item #8206, 12/package

**Reusable Callus Cushions**

- Item #8207, 6/package

**Metatarsal Support Pads (See page 37)**

- Item #8208, 1 pair/package

**Bunion Shield**

- Item #8212, 1/package

*"I have to say the Pedi-Gel pads have been a big hit with the patients and a really high quality product. Truly amazing how long they last. The shapes I find most helpful are the U-shape Pad to off-load painful IPK, warts, pre-ulcerative lesions, and the Met Pad to treat general metatarsalgia and lesser MTPJ pain."*

*Dr. W., South Hadley, MA*



Ball-of-Foot Pads

Arch Pads

Reusable Callus Cushions



Callus Pads

Dancer's Pads

Corn Cushions



Bunion Shield



Sheets (2" x 3 1/2")

**New!**

# FELTastic® Felt Foot Pads

## Are all Felts the same?

No. Ours is FELTastic! Higher cotton and wool content feels better and lasts longer. The proper mix of Rayon adds durability. Compare our medical-grade FELTastic® Felt to others — feel the difference! Your customers will, too.

All these self-adhesive pads will stick to feet, in shoes, on, or under insoles to relieve pain.



#2105



#2105

#8504

## **New! FELTastic® Heel Pads** *Protect & Lift to Relieve Pain*

These absorb impact shock, raise the heel to relieve stress on the plantar fascia and Achilles, and compensate for leg-length discrepancies that cause imbalance in hips and back. 100 per pack. **#8534**



#8534

## **New! FELTastic® Arch Pads** *Instant Support Makes Shoes & Feet Feel Better*

Relieve tired, achy and flat feet, limit pronation and add extra support to footwear lacking an arch raise. Help ease heel, arch and ball-of-foot pain by redistributing weight and pressure. One scissor-trimmable size fits most. 10 per pack. **#8536**



P8536

## **New! FELTastic® Flat Metatarsal Pads**

*Stick to Shoes or Feet for Immediate Comfort*

Instant relief for ball-of-foot pain, calluses and sore spots. Cushion, protect and offload sensitive areas. One scissor-trimmable size fits most. 100 per pack.

**#8501** 1/8" thickness

**#8502** 1/4" thickness



#8501/#8502

## **New! FELTastic® Contoured Metatarsal Pads** *Shaped to Support & Protect*

These ball-of-foot pads feature a rounded, 'skived' shape that supports the metatarsal arch and bones to improve foot structure and biomechanics. Relieve ball-of-foot discomfort, calluses, pressure and sore spots, even stubborn neuroma pain, and more. Apply to feet, in shoes, atop or under insoles. Sizing is approximate.

Small Fits Women's shoe sizes 4-8, Men's 5-9. Large fits Women's shoe sizes 9+, Men's 10+.

Small **#8530-S**

Large **#8530-L**



#8530-S



#8530-L

## **New! FELTastic® Oval Callus Pads** *Surround & Protect Tender Areas*


Oval shape helps offload calluses and sensitive spots to relieve pressure and pain. Frequently applied to ball-of-foot, also heels. Use 'as is' or cut in half and place on each side of a larger sore spot. One size. 100 per pack.

**#8504**



#8504





**Description**

FELTastic® Self-Adhesive Corn Pads, Beige, 1/8"

FELTastic® Self-Adhesive Callus Pads, U-Shaped, Beige, 1/8"

FELTastic® Self-Adhesive Callus Pads, Oval Shaped, Beige, 1/8"

Pedi-Patch™ Self-Adhesive Moleskin Foot Protection Pads

FELTastic® Self-Adhesive Flat Metatarsal Pads, Beige, 1/8"

FELTastic® Self-Adhesive Flat Metatarsal Pads, Beige, 1/4"

FELTastic® Self-Adhesive Contoured Metatarsal Pads, White

FELTastic® Self-Adhesive Contoured Metatarsal Pads, White

FELTastic® Self-Adhesive Flat Heel Pads, White, 1/4"

FELTastic® Self-Adhesive Flat Heel Pads, White, 1/4"

FELTastic® Self-Adhesive Flat Heel Pads, White, 1/4"

FELTastic® Self-Adhesive Arch Pads, White

Self-Adhesive Felt Roll, White, 1/8" x 6" x 2.5 yd.

Self-Adhesive Felt Roll, White, 1/4" x 6" x 2.5 yd.

Self-Adhesive Foam Roll, White, 1/8" x 6" x 2 yd.

Self-Adhesive Foam Roll, White, 1/4" x 6" x 2 yd.

**Size**

One Size

One Size

One Size

One Size

One Size

One Size

Small

Large

Medium

Large

XLarge

One Size

**Product#/Pack Bulk**

#8005-100

#8023-100

#8504-100

#2105-100

#8501-100

#8502-100

#8530-S-10

#8530-L-10

#8534-M-100

#8534-L-100

#8534-XL-100

#8536-10

#8520

#8521

#8620

#8621

**OTC#/Pack**

#P15/10 pack

#P16/8 pack

#P8504/8 pack

#P2105/6 pack

#P8501/2 pack

#P8502/2 pack

#P8530-S/2 pack

#P8530-L/2 pack

#P8534-M/2 pack

#P8534-L/2 pack

#P8534-XL/2 pack

#P8536/2 pack

## PediPlast® Custom Digital Devices in Less Than 5 Minutes!



*Just a few of the devices you can make in minutes — for accommodation of digital problems or post-surgical alignment and compression.*

### **Clean, fast and easy**

PediPlast® is a putty-like material that you mold directly to your patient's foot, where it cures to a soft, permanent form. You control the density and firmness of the finished devices by adjusting the amount of hardening agent added. The process takes only minutes, is clean and trouble-free. After just a few uses, you'll be confident with the material, the process and the consistency of the resulting appliances.

***Keep PediPlast® handy in your office for both routine and difficult cases not easily treated with other methods.***

### **Soft, yet durable**

Patients appreciate PediPlast® footcare appliances because they're custom made just for them, and the devices feel good and look natural. The material is non-toxic, non-irritating, hypo-allergenic, soft, yet durable and tear resistant. Washable with soap and water, too. No wasted material, no smell or residue, no refrigeration, no special handling required. PediPlast® has an indefinite shelf life and will not harden until you add the curing agent. Scissor-trim or grind finished devices to make additional alterations after hardening, if needed.

The possibilities are truly unlimited, giving you the powerful ability to respond to situations as they arise — with a quick, custom-made solution to most digital problems — either temporary or permanent.

PediPlast® is a moldable silicone compound that lets you fabricate soft, durable, custom digital devices in less than 5 minutes. It's fast, easy and versatile, providing you with the ability to relieve a variety of digital conditions.

Use PediPlast® to:


- Cushion and Protect
- Align and Separate

Use PediPlast® to treat dozens of foot conditions, including:

- Bunions and Corns
- Hammer, Mallet, Claw, Crooked, Overlapping Toes
- Amputations
- And More!

### **Exceptional value, too.**

PediPlast® is economical. Your cost is usually less than \$2.00 per device, and the 250 gram jar can produce as many as 50 appliances, depending on their size and complexity. The cost is affordable for most patients — typically from \$10 to \$25 per appliance.



*Custom Hammer Toe Crest, Custom Toe Separator, Just mix & mold!*



**Just Mix, Mold & Apply!**



800g



400g



250g



100g

**Description**

- 100 grams jar & 2 tubes of curing agent
- 250 grams jar & 3 tubes of curing agent
- 400 grams jar & 4 tubes of curing agent
- 800 grams jar & 10 tubes of curing agent

**Item#**

- #8000**
- #8001**
- #8002**
- #8003**



#8000



#8001



#8002



#8003



**Directions for use:**

1. Combine PediPlast with 1cm of curing agent per one scoop of compound. Add slightly less curing agent for softer devices that cure slowly. Add slightly more curing agent for firmer devices that cure faster.



2. Mix thoroughly in hand for approximately 20 seconds. (Notice changing density.)  
 3. Roll into ball, then into an elongated shape before molding to patient's foot. (This helps to avoid seams in the final appliance.)



4. Apply to patient's foot and mold until desired shape is achieved.  
 5. In 4-6 minutes, PediPlast will reach its permanent form. Alterations can be made — cut or grind after permanent form is reached.

**Don't be fooled by imitations, insist on PediFix® quality products**

**Hammer Toe Crests**

Relieve Toe Tip Pressure, Support Bent-Under Toes, Make Walking More Comfortable

- Relieve toe tip corns and irritations
- Cushion and support hammer, claw, mallet, arthritic and bent-under toes
- Ease metatarsal and forefoot pain
- Anatomically-correct shape for maximum comfort and effectiveness

**Choose Suede, Felt or Visco-GEL®**

Suede and felt-covered Toe Crests feature soft, durable fabrics, and adjustable elastic bands that hold them in place. Visco-GEL® Crests feature our proprietary Gel material and a flexible, fixed loop.



Suede & Felt Crest toe loops are adjustable. Gel Toe loops are flexible, but fixed.



3/pk	Deluxe Suede	Felt	Visco-GEL®
Small Left	#8154A-SL	#8154B-SL	#3026
Small Right	#8154A-SR	#8154B-SR	#3025
Medium Left	#8154A-ML	#8154B-ML	#3028
Medium Right	#8154A-MR	#8154B-MR	#3027
Large Left	#8154A-LL	#8154B-LL	#3030
Large Right	#8154A-LR	#8154B-LR	#3029
Extra Large Left	#8154A-XLL	#8154B-XLL	
Extra Large Right	#8154A-XLR	#8154B-XLR	
Universal-OSFM	N/A	N/A	#3035

**Hammer Toe Crests Sizing Chart**

Small fits Women's Shoe Sizes 4-7  
 Medium fits Women's Shoe Sizes 8-10, Men's 6-8  
 Large fits Women's Shoe Sizes 11+, Men's 9-10  
**New!** Extra Large fits Men's Shoe Sizes 11+

**1/pk Item#**

Small Left	#1036-SL
Small Right	#1036-SR
Large Left	#1036-LL
Large Right	#1036-LR

*"Your Hammer Toe Crests are an excellent method for relieving pain associated with corns and arthritis related to hammer toe deformities"*

Dr. R.R. Jupiter, MD

*"I continually save valuable time by not having to make palliative products. I get excellent results using the PEDIFIX product line."*

Dr. L., San Carlos, CA

**New!**

**Visco-GEL® B2Splint™**

*Patented Toe Spacer/Crest Combination*

This new device aligns & separates 1st and 2nd toes to relieve Hallux Valgus irritation. Gel loop limits Crossover Toe, and the Sulcus support Toe Crest lifts hammer, claw and mallet toes and tips. Effective pre-op, post-op and non-operatively. Test this before surgery. Choose Right or Left, Small fits shoe sizes up to 8, Large fits 9+. OTC Name: ToeCoach™

1/pk	Item#
Small Right	#2853-SR
Small Left	#2853-SL
Large Right	#2853-LR
Large Left	#2853-LL



Spacer separates 1st and 2nd Toes

Loop limits under or overlapping, 'crossover toe', keeps device in place



Sulcus extension supports Hammer Toes, Offloads Toe Tips

**Visco-GEL® Hammer ToeCrutch®**

*Relief for Hammer or 'Claw' Toes Innovative New Single-Toe Design*

Specially contoured to support the sulcus, elevate toe tip and relieve distal pressure. Helps involve hammer and claw toes in gait for improved biomechanics. Soft, extra-wide, Visco-Gel® toe loop protects, hydrates & softens typical 'PIP' toe-top corns, and keeps device in place.

4/pk	Item#
Small	#1037-S
Medium	#1037-M
Large	#1037-L



**Hammer ToeCrutch® Sizing Chart**

Small fits smaller toe (small buttress, 13 mm diameter Gel toe loop, 28 mm length)  
 Medium fits average toes (medium buttress, 14 mm diameter Gel toe loop, 32mm length)  
 Large fits big toes (large buttress, 16 mm diameter Gel toe loop, 37 mm length)



### Budin Toe Straighteners

(Budin Toe Splints)

Toe Straighteners Align Problem Toes — Conservatively or Post-op

- Straighten crooked, overlapping digits
- Encourage hammer toes to ‘flex’ to help postpone rigidity
- Soft elastic loops guide deformed toes into place
- Ball-of-foot cushion helps protect metatarsal heads
- Plastic disk (hidden inside) stabilizes loops
- Effective post-operative splint to speed healing

*Ideal after Capsulotomy or Tenotomy procedures*



1/pk	Item#
Single Toe - OSFM*	<b>#8155</b>
Double Toe - OSFM*	<b>#8157</b>
Triple Toe - Left	<b>#8158-L</b>
Triple Toe - Right	<b>#8158-R</b>

\*interchangeable for right or left foot

This product contains natural rubber latex which *may* cause allergic reactions.




### Comfort ToeWraps™

*Align Crooked, Overlapping, Injured or ‘Hammer’ Toes*

*Wrap-Around Cushioning Relieves Rubbing & Irritation*

- Adjustable Loops Separate & ‘Buddy’ Adjacent Toes
- Relieves Between-Toe Pain
- Gently Stabilize & Reposition Without Irritation
- Comfortable, Convenient & Long-Lasting
- Effective as a Toe Splint for Injured, Broken or Post-Op Recovery



6/pk	Item#
OSFM	<b>#6055</b>

Universal for Left or Right Foot

### Velcro® Digit Wraps

Soft, Pre-Made ‘I’ Splint for Adjacent Toes

- Align crooked, overlapping digits
- Straighten flexible hammer toes
- Effective conservatively, post-surgically or following injury
- Easily removed and reapplied for bathing, therapy or inspection

6/pk	Item#
OSFM	<b>#6051</b>



### ToeSnug® Adjustable Hammer Toe Wrap™

*Relieves Hammer Toe Pain*

Versatile, new device, cleverly engineered for multiple uses. Supports, separates, aligns and comforts crooked, overlapping and claw toes, cushions hard and soft corns, relieves Metatarsalgia, too. Customize to fit around individual or groups of digits. Dual-action design takes pressure off toe tips and encourages flexion and extension of hammer toes. Patent pending. Interchangeable for right or left foot. Sizing is approximate. Small fits shoe sizes 4-7. Medium fits shoe sizes 7-10. Large fits shoe sizes 10+.



Before



After

1/pk	Item#
Small	<b>#8163-S</b>
Medium	<b>#8163-M</b>
Large	<b>#8163-L</b>




**Visco-GEL® All-Gel Digital Caps**

Thin side wall design takes up less space in shoes, 'Comfort Vent' for breathability, extra distal padding for added toe tip protection.

Size	10 Pack Item#
Small/Medium (15 mm diam.)	#1075
Large/XLarge (20 mm diam.)	#1076
XXLarge (25 mm diam.)	#1077



Stretches easily



*"...Your Visco-Gel products are really great — consistently good quality, expert designs, and my patients love them. Your customer service people are great to deal with — always helpful and very knowledgeable..."*

*Dr. H., Katonah, NY*



**New!**

**Visco-GEL® Heavy Duty 'Sports' Toe Cap**

*Heavy-Duty, Last Longer, Resist Tearing, Easier for Elderly Patients to Find*

10/pk	Item#
Small (1/2" diam.)	#1078-S
Large (3/4" diam.)	#1078-L



Hand Sewn Durable Cover

**Nylon-Covered Foam Toe Caps**

*Instead of Gel*

- Proprietary foam
- Extra distal padding for added toe-tip protection
- Our exclusive, hand-sewn nylon cover adds durability and reduces friction


3/pk	Item#
Small (1/2" diam.)	#8134-S
Medium (3/4" diam.)	#8134-M
Large (1" diam.)	#8134-L

**Visco-GEL®**

**Fabric-Covered Digital Caps**

Choose fabric sleeve: 'Ribbed' is most compressive. 'Mesh' offers more stretch. 'Knit' style features no seam, and is least compressive.

Ribbed	Seamless Knit	Mesh	6/pk
#1050	#1053	#1056	Small/Medium (15 mm diam.)
#1051	#1054	#1057	Large/XLarge (20 mm diam.)
#1052	#1055	N/A	XXLarge (25 mm diam.)



#1050  
#1053  
#1056  
S/M

#1051  
#1054  
#1057  
L/XL

#1052  
#1055  
XXL



Ribbed



Seamless Knit



Mesh



Kit-2



Kit-1

**Visco-GEL® "Immediate Relief" Treatment Kits**

**29 & 62 Piece Assortments**

Every office should have one! Compact, portable kits keep a variety of Visco-GEL® Foot Protection Products handy and organized. You'll treat dozens of common foot injuries, ailments and deformities, and provide patients with instant pain relief. Re-stock anytime. Return unused product for exchange, refund or credit.

**#Kit-1  
Standard Kit Contents**

- 4 Corn Pads
- 3 Ribbed Digital Caps
- 2 Mesh Digital Caps
- 3 All-Gel Digital Caps
- 6 Toe Separators
- 3 Toe Spacers
- 2 Fully-Coated Ribbed Tubes
- 2 Fully-Coated Mesh tubes
- 2 Bunion Guards
- 2 Universal Toe Crests

**10 Products, 29 Pieces**

**#Kit-2  
Deluxe Kit Contents**

- 6 Corn Pads
- 4 Toe Spacers with Loop
- 6 Ribbed Digital Caps
- 6 All-Gel Toe Caps
- 9 Toe Separators
- 9 Toe Spacers
- 2 Fully-Coated Ribbed Tubes
- 2 Fully-Coated Mesh Tubes
- 2 Toe Spacer/Bunion Guards
- 4 Bunion Guards
- 2 Universal Toe Crests
- 2 Gel Sheet w/Cover
- 2 Gel Ball-of-Foot Pads
- 2 Gel Arch Pads
- 4 Gel Callus Pads

**15 Products, 62 Pieces**

*Kit contents and prices are subject to change.*

## Digit Care Made Easy



Soft, Durable  
3-Layer Foam

Large  
Medium  
Small


*"After wearing them for two months, to my delight, the ulcer healed and the Toe Separators saved an amputation of my right toe"*

Mrs. L.B. Windsor, Ontario

### 3-Layer Toe Separators™

Soft, soothing, exclusive PolyFoam™ outer layers, slightly firmer middle layer.

12/pk	Item#
Small	#8130-S
Medium	#8130-M
Large	#8130-L
XLarge	#8130-X
6/pk-MIX	Item#
2S, 2M, 2L	#8130-MIX/6




Size	Length	Thickness
Small	30mm	10mm
Medium	42mm	12mm
Large	50mm	13mm
XLarge	50mm	16mm

**New!**

### Moleskin-Lined Toe Separators

Soft, soothing, moisture-absorbing Moleskin outer layer, foam middle layer. Comfortable. Durable. Washable.



12/pk	Item#
Small	#8131-S
Medium	#8131-M
Large	#8131-L
XLarge	#8131-X

### Tubular Foam Toe Bandages™

Soft, fabric-lined foam, cut to size

6/pk (12" long ea.)	Item#
Small (5/8" diam.)	#8137-5
Medium (3/4" diam.)	#8137-3
Large (1" diam.)	#8137-1



### Visco-GEL® Toe Separators

Proprietary Gel divides toes gently, soothes and softens.

15/pk	Item#
Small	#1100
Large	#1102
XLarge	#1103



### PolyFoam™ Toe 'Comb' Cushion

- Super soft foam in a unique anatomical shape.
- Use as-is or trim for use between 2 or 3 digits.
- Position above or underneath toes, to relieve pressure on toe tops or tips

12/pk	Item#
OSFA	#8230



### Visco-GEL® 5-ToeBuddy®

Separate, Align, Stretch Toes

These multi-toe spacers relieve toe-to-toe irritations, crooked, overlapping digits, between-toe corns, and more. Gently align and separate, inside shoes and socks, or barefoot for a relaxing, all toes stretch. Anatomically-shaped to fit Right and Left feet precisely. Scissor-trimmable for easy customization for fewer toes. 1 pair.

1 pr/pk	Item#
	#P1114



1st Toe Loop keeps device in position



Contoured Spacers gently separate crooked toes



**Visco-GEL® Foot Protection**



**Visco-GEL® Corn Pads**

Soft Ribbed fabric, **Half-Coated** with Gel

12/pk	Item#
Small/Medium	<b>#1000</b>
Large/XLarge	<b>#1001</b>

**New!**

**Visco-GEL® Rolled Digital Corn Pad**

Soft, ribbed fabric with Gel dot.



12/pk	Item#
Mini	<b>#1010-MINI</b>
Small/Medium	<b>#1010-M/L</b>
Large/XLarge	<b>#1010-L/XL</b>



**Visco-GEL® Digital Tubes**

Thin wall design fits comfortably in shoes.

10/pk	Item#
Small (15 mm diam.)	<b>#1027</b>
Medium (20 mm diam.)	<b>#1028</b>
Large (25 mm diam.)	<b>#1029</b>



**Visco-GEL® Fully-Coated Digital Tubes**

Soft Mesh or Ribbed fabric, with Gel lining. Cut to size

3/pk (6" long ea)	Mesh	Ribbed
Narrow (15 mm diam.)	<b>#1202</b>	<b>#1200</b>
Wide (20 mm diam.)	<b>#1207</b>	<b>#1201</b>



**Visco-GEL® Little Toe Sleeves™**

Specialized size for smaller toes. Conical, tapered shape won't roll or migrate.

10/pk	Item#
OSFM	<b>#1026</b>



**Visco-GEL® Half-Coated Digital Strips**

Soft Knit fabric, with Gel lining, cut to size

2/pk (24" long ea)	Item#
Small/Medium (15 mm diam.)	<b>#1150</b>
Large/XLarge (20 mm diam.)	<b>#1151</b>



**Visco-GEL® Digital Pad Strip**

Soft Knit fabric, with 1" diameter Gel dots, cut to size

2/pk (24" long ea)	Item#
Mini (10 mm diam.)	<b>#1177</b>
Small/Medium (15 mm diam.)	<b>#1175</b>
Large/XLarge (20 mm diam.)	<b>#1176</b>

## Help for Crooked, Overlapping Toes



### Visco-GEL® ToeBuddy®

Two soft Gel toe loops keep this Spacer in perfect position between digits to separate and help align crooked, overlapping or otherwise deformed first and second toes. One size fits most.

**4/pk**    **Item#**  
OSFM    **#1033**



### Visco-GEL® Little ToeBuddy®

*Separates & Aligns Crooked and Overlapping Toes*

Contoured Spacer fits comfortably between smaller toes to absorb pressure and friction and relieve pain. Cushions corns, divides and aligns crooked and overlapping toes. Soft, flexible loops slip over toes to keep Spacer in position. Visco-GEL® soothes and softens skin. One size fits most.



*Works on 4th and 5th or 3rd and 4th toes*



*New Wider Toe Loops*

**4/pk**    **Item#**  
OSFM    **#1035**



**New!**

### Visco-GEL® ToeBuddy® 'Long'

*Ideal for longer toes & more severe 1st and 2nd toe deformities*

Our new design features a longer, extended Toe Spacer to more effectively separate digits. Ideal for larger-sized feet, and more severe 1st and 2nd toe deformities.

**4/pk**    **Item#**  
OSFM    **#1034**



*Hallux Limitus Application*



### Visco-GEL® Toe Spacer with Fabric Sleeve

Separates toes, keeps Spacer in its proper position. Fits on first or second toes. Reversible, can be applied with Gel 'in' or 'out'

**4/pk**    **Item#**  
Medium    **#1006-M**  
Large    **#1006-L**





**Treatment Tip:**  
Position with the Flex Slit facing down, which allows these to deform efficiently on impact.



**Visco-GEL®  
Stay-Put Toe Spacers®**

*Gel Loop Keeps Spacer in Position Between 1st & 2nd Toes*

Soft Gel loop fits around first or second digit to keep Spacer in position between toes.

4/pk	Item#
Medium	#1031-M
Large	#1031-L



Size	Height
Medium	28mm
Large	30mm

**Visco-GEL® Toe Spacers**

Mild flanges help keep these soft Gel spacers in place between 1st & 2nd toes.

10/pk	Item#
Small	#1125-10
Medium	#1126-10
Large	#1127-10

*"I have a bad case of bunions but can't have surgery... these Visco-GEL® Toe Spacers are a life saver... I am so thankful for them... this has relieved so much of my pain, with these I can stay comfortable. Thanks for a wonderful product"* Victoria from Vermont



**Original Toe Spreaders**

Firmer than Gel, with a more pronounced flange to better separate and align toes, conservatively or post-op.

This product contains natural rubber latex which may cause allergic reactions.

12/pk	Item#
Small	#8129-S
Medium	#8129-M
Large	#8129-L

**New!**

**Visco-GEL® Stay-Put Toe Separators® for Smaller Toes**

*Gel Loop Keeps Spacer in Position Between Toes*

Soft contoured Gel Separators absorb pressure and friction to relieve between-toe pain. Soothe corns and gently divide toes that rub. Help align and separate crooked, overlapping digits. 'Hourglass' shape for all day comfort. Soft Gel loop fits around any small toe to keep Separator in position.



10/pk	Item#
Small	#1104-S
Medium	#1104-M
Large	#1104-L

Size	Height
Small	20mm
Medium	25mm
Large	29mm

## Innovation in Personalized Orthodigital Care



#1030

**New!**

### Visco-GEL® Customizables®

*Adjustable Relief, Comfort & Shoe Fit  
Make You The Hero*

Now you can adjust the effect of our Gel devices to each patient's needs and shoe fit constraints. Patent-pending Customizables™ allow personalized relief for common toe ailments and deformities. Available only to healthcare professionals and their patients, so you can guide and adjust the fit and function of each device.

- Choose the appropriate device for your patient's condition
- Personalize with No, Soft or Firm durometer PediPlug™ insert. Test for comfort and shoe fit.
- Show Patients how to exchange Pedi-Plugs™ by stretching the device opening, removing and inserting different density inserts.
- Encourage Patients to gradually train the foot to accept firmer durometer Plugs.

Enjoy the benefits of customized comfort!

*Customizables® make our original designs even better!*

### Customizables® Stay-Put Toe Spreader®

**Condition:** 1st and 2nd Toes  
Interdigital Irritation, Overlapping, Bunions.  
Gently divides, separates toes that rub. Aligns big toe to reduce pressure on the 1st MPJ prominence. Soft, flexible loop slips over first or second toes to keep Spreader in perfect position. Interchangeable for right or left foot. Choose firm, soft or no Pedi-Plug depending on ailment severity and shoe fit.

1/pk

Small/Medium

#1030-S/M-CTZ

Large/XLarge

#1030-L/XL-CTZ



### Customizables® Toe Spacer/Bunion Guard Combo

**Condition:** Bunions & Overlapping Big Toe  
Contoured Gel shield protects 1st MPJ bunion prominence. Attached Separator divides, aligns big and second toes. Choose Firm, Soft or No PediPlug™ Insert for personalized comfort and shoe fit.

1/pk

OSFM

#1315-CTZ



## Customizables® Hammer ToeCrutch®

**Condition:** Hammer, Claw, Mallet, Contracted Toes

Elevate toe-tips to relieve distal pressure. Helps involve bent-under toes in gait for improved biomechanics. Soft, extra-wide Visco-GEL® toe loop protects, hydrates and softens typical toe-top corns and keeps device in place. Personalize toe lift with firm, soft or no PediPlug™ insert.

- 1/pk**
- Small #1037-S-CTZ
  - Medium #1037-M-CTZ
  - Large #1037-L-CTZ



## 10 Seconds to a Customized Fit



**Option A: Soft Fit**  
Use no insert for the most gentle fit & a thinner profile in footwear



**Option B: Medium Fit**  
Use the medium insert for a softer separation between toes



**Option C: Firm Fit**  
Use the firm insert for a stronger separation between toes

## Customizables® ToeBuddy®

**Condition:** Crooked, Overlapping Toes, Bunions, Between-Toe Corns.

Two soft Gel loops position anatomically-shaped spacers between toes to absorb pressure and friction, divide and align. Big Toe size fits between first and second toes, helps position the big toe to reduce bunion prominence. Little Toe size fits between lesser digits, most commonly 2nd and 3rd, or 4th and 5th. Personalize the fit with Firm, Soft or no Pedi-Plugs™.

- 1/pk**
- Little Toe #1035-CTZ
  - Big Toe #1034-CTZ



## AdjustaLoop™ Toe Crest

**Condition:** Rigid Hammer, Claw, Mallet, Contracted Toes

This exclusive Visco-GEL® Toe Crest is positioned with a soft, elasticized band that is uniquely adjustable on both its ends to help prevent patient overtightening. Supporting the sulcus, the device relieves pressure on toe tips and distal irritations, and involves the toes in ambulation to make walking easier for patients. Proper anatomical curvature for comfort and effectiveness. Choose Right or Left, S for shoe sizes approximately 4-8, L for sizes 9+.

- 1/pk**
- Small Right #3037-SR
  - Large Left #3037-SL
  - Small Right #3037-LR
  - Large Left #3037-LL



[www.pedifix.com/t-Customizables.aspx](http://www.pedifix.com/t-Customizables.aspx)

Pre-op, Post-op, or Non-Surgical Relief



Contoured shape fits more comfortably



Ultra-thin Design



Stretchable Toe Loop

**Visco-GEL® Bunion Guards**

Relieves Pain, Cushions & Protects Sensitive MP Joint

- Absorbs shoe pressure, friction
- Proprietary Gel soothes and softens
- Thin profile design fits comfortably in most footwear
- Alternate positioning — under sesamoids



“...Ideal protective device for bony prominence, with thin profile...”

Dr. Nicholas Romansky, DPM

1/pk	Item#
Hallux - OSFM	#1316
Tailor's - OSFM	#1319



**Visco-GEL® Dual-Action Toe Spacer/Bunion Guard Combo**

Contoured shape fits more comfortably

Cushions and protects 1st metatarsal prominence and prevents irritation between first and second toes.

1/pk	Item#
OSFM	#1315



**New!**

**Visco-GEL® ToeBuddy® Bunion Guard**

Loop Helps 2nd Toe Alignment

1038

Now with SmartGel® Technology!

Dual-action device aligns first & second toes, protects bunion joint. Contoured Gel cup cushions the first metatarsal-phalangeal head from shoe pressure and friction, while dual toe loops help separate and align first and second toes. New, improved design. Comfortable. Durable. Washable.

1/pk	Item#
OSFM	#1038



**Visco-GEL® Forefoot Protection Sleeve**

Wraps Feet in Gel Comfort

Protective Visco-GEL® cushioning extends from the bunion joint, across the ball-of-foot and around the 5th toe to keep sensitive feet comfortable. Absorbs pressure to relieve bunion and metatarsal pain. Reversible to fit right or left foot. Washable, reusable, effective for months.

Soft'N Stretch™ Fabric

1/pk	Item#
Small/Medium	#1455-S/M
Large/XLarge	#1455-L/XL



Hallux Protection Tailor's Protection



**New!**

**Visco-GEL® Little ToeBuddy® Bunion Guard**

1/pk	Item#
OSFM	#1320

1320



**Visco-GEL®  
Bunion Care™ Relief Sleeve**

Cushion, Protect, Absorb Pressure & Friction

- Relieve pain, makes shoes comfortable
- Soft'N Stretch™ fabric cover positions Visco-GEL® pad over MP joint
- 'Thin Dress' model features 3mm thick Gel pad and fits in most fashion footwear
- 'Extra Protection' model features 5mm thick pad that requires more space in shoes
- Uncovered style keeps Gel in contact with skin surface to soothe, soften
- Covered version places fabric next to skin
- **New Bunion Care™ Relief Sleeve with Aperture** — protects more severe bunion deformities

1/pk (3mm pad) "Thin Dress"	Covered	Uncovered
Small/Medium	#1303	#1306
Large/XLarge	#1304	#1307
XXLarge	#1308	
1/pk (5mm pad) "Extra Protection"		
Small/Medium	#1309	
Large/XLarge	#1310	#1313
1/pk With Aperture		
Small		#1311-S
Large		#1311-L



**FELTastic® Bunion Shield**

Soothing MP Joint Protection Fits Comfortably In Most Shoes

- 2 layers of soft, durable felt reduce pressure, friction.
- Elastic loop keeps Shield in position

1/pk	Item#
OSFM	#8148



**Visco-GEL®  
Forefoot Protection Sockette**

Complete plantar forefoot protection, from toe tips to metatarsal heads. Thin layer of Gel covered on both sides with fabric.

1/pk	Item#
Small/Medium	#1342
Large/XLarge	#1343

interchangeable for right or left foot



**Pedi-GEL® Bunion Shield**

Self-adheres to skin to cushion, protect and soothe. Washable and reusable.

1/pk	Item#
OSFM	#8212



## FootSplints™ for Post-Surgical Management



**New!**

### Comfort Bunion Splint™

*Aligns Toe, Protects 1st MTPJ*

Soft, elasticized fabrics gently align and position the big toe for improved comfort or post-op alignment. Velcro® attachments allow complete customization. Gel pad cushions the 1st metatarsal head prominence. Low friction surfaces minimize shear forces.

1/pk	Item#
Small Left	#6025-SL
Small Right	#6025-SR
Large Left	#6025-LL
Large Right	#6025-LR

**New!**

### HV Hinge™ Bunion Splint

*Allows Flexion/Extension Motion, Controls Alignment*

While it straightens the big toe, the rotating axis of this unique splint allows plantar and dorsal motion, which many surgeons find advantageous at select times during the post-op recovery process. Simple, and low cost, they can play an effective role in transitioning your patients to increased ranges of motion and mobility, even if used only for a short time. Intended for non-ambulatory use. One size fits most, fits either foot, flip for Right/Left use.



1/pk	Item#
OSFM	#6038

### Hallux Valgus SoftSplint™

Comfortable, Effective Post-Op Bunion Splinting

- For ambulatory and non-ambulatory use
- Maintains ideal hallux positioning and constant MP joint alignment
- Completely adjustable positioning control
- Soft fabric construction

1/pk	Item#
Small Left	#6026-SL
Small Right	#6026-SR
Medium Left	#6026-ML
Medium Right	#6026-MR
Large Left	#6026-LL
Large Right	#6026-LR

HCPCS Code: L3100

*"I use the Nighttime Bunion Regulator in patients who want to postpone surgery with good results. Post-operatively, they eliminate time spent making and changing homemade splints and the inconvenience they cause my patients"*

J.W. MD

*"I find using the NightSplint when bunions just begin, it avoids the need for surgery. Really very convenient"*

J.P. MD

### FootSplints Sizing Chart

Small generally fits Shoe Sizes 4-7  
 Medium generally fits Shoe Sizes 8-10  
 Large generally fits Shoe Sizes 10+  
 Exact sizing is not critical for these products.

### Hallux Valgus NightSplint™

Non-Ambulatory Positioning Control

- Immobilizes hallux to speed healing
- Adjustable varus-valgus alignment to maintain ideal surgical result

1/pk	Item#
Small Left	#6035-SL
Small Right	#6035-SR
Medium Left	#6035-ML
Medium Right	#6035-MR
Large Left	#6035-LL
Large Right	#6035-LR

*Also comforting conservatively for patients who don't want surgery*





**Forefoot Compression Sleeve™**

- Controls and Reduces Edema After Injury or Surgery
- Provides Grade II Compression (20-30 mmHg)

1/pk	Item#
Small (fits women's 6-8, men's 4-6)	#6027-S
Medium (fits women's 9-11, men's 7-9)	#6027-M
Large (fits women's 12+, men's 10-12)	#6027-L

See our Compression Ankle and additional edema reduction products on page 33.



**Reduces Post-Op Swelling**



**Darco® Flexible Post-Op Shoe**

Reduces Stress on Surgical & Fracture Sites to Speed Healing

- Semi-flexible sole conforms to foot contours throughout the gait cycle, makes ambulation easier and more comfortable.
- Soft nylon upper, self-molding insoles protect foot surface.
- Interchangeable right or left foot.

1/pk	Item#
Women's Small (fits 4-6)	#6210-L-S
Women's Medium (fits 6½-8)	#6210-L-M
Women's Large (fits 8½-10)	#6210-L-L
Men's Small (fits 6-8)	#6210-M-S
Men's Medium (fits 8½-10)	#6210-M-M
Men's Large (fits 10½-12)	#6210-M-L
Men's XLarge (fits 12½-14)	#6210-M-XL
Children's (fits 12-1)	#6210-Y

Order 12-49 for just \$9.00 ea.  
 Order 50+ for just \$7.70 ea.  
 HCPCS Code: L3260



**Hallux Valgus DaySplint™**


Adjustable Alignment Control

- Simple varus-valgus adjustment
- Permits hallux flexion-extension

1/pk	Item#
Small Left	#6022A-SL
Small Right	#6022A-SR
Large Left	#6022A-LL
Large Right	#6022A-LR

Small fits narrow feet. Large fits wide feet.

## Visco-GEL® Forefoot Protection



Most popular version covered inside and out

*"I went to see my doctor and he said try this Universal Metatarsal Strap...it is a wonderful product... thank you very much!!!"*

*Louise from Hawaii*

Also available in uncovered version with moisturizing Visco-GEL® against the skin.



### Visco-GEL® Universal Met Strap

Soft'N Stretch™ fabric forefoot sleeve with Gel pad inside. Protects met heads from pressure.

1/pk	Covered	Uncovered
Small/Medium Right	#1225	#1229
Small/Medium Left	#1226	#1230
Large/XLarge Right	#1227	#1231
Large/XLarge Left	#1228	#1232



Thin design, comfortable in everyday footwear

### Visco-GEL® Silicone Thin Forefoot Cushion

2mm thick Gel with Soft'N Stretch fabric forefoot wrap and 1st toe loop.

1/pk	Item#
Small Right	#1321
Small Left	#1322
Large Right	#1323
Large Left	#1324



### Visco-GEL® Metatarsal Wrap

*Helps Relieve Ball of Foot Pain*

Stretchable fabric positions a soft silicone Gel cushion that protects metatarsal heads. For use under or over socks, or barefoot around the house. Thin profile fits comfortably in most footwear. Washable and reusable for months. Interchangeable – fits right or left foot.


1/pk	Item#
Small/Medium (fits narrow or average width feet)	#1295-S
Large/XLarge (fits wide feet)	#1295-L

### Arch Support Bandages

Compressive elastic midfoot bandage gently stabilizes, repositions metatarsal shafts, heads

- Supports plantar fascia to help relieve heel pain
- Eases ball of foot pain, bunion and Tailor's bunion pain
- Approximately 50 mm Hg compression

2/pk	Item#
Small (fits midfoot circumference 18-19 cm)	#6000-19
Medium (fits midfoot circumference 20-21 cm)	#6000-21
Large (fits midfoot circumference 22-23 cm)	#6000-23
XLarge (fits midfoot circumference 24-25 cm)	#6000-25



### Arch Binders™ with Metatarsal Pad

Compressive elastic midfoot bandage with removable metatarsal pad repositions, stabilizes and supports metatarsal heads, shafts.

2/pk	Item#
Small (fits midfoot circumference 18-19 cm)	#6002-19
Medium (fits midfoot circumference 20-21 cm)	#6002-21
Large (fits midfoot circumference 22-23 cm)	#6002-23

interchangeable for right or left foot



### Sesamoid Relief Sleeve™ Improved Design!

Lightweight Lycra® forefoot binder positions a soft foam ball-of-foot pad with sesamoid relief aperture under the first metatarsal head to limit pressure and shock and relieve pain. Washable and reusable. Sizing is approximate. Small fits narrow and average width feet. Large fits wide feet.

1/pk	Item#
Small Left	#6003-SL
Small Right	#6003-SR
Large Left	#6003-LL
Large Right	#6003-LR



Covered gel pad

### Super Metatarsal Cushion

Thick PolyFoam™ padding on soft elastic midfoot sleeve.

1/pk	Item#
Small (fits midfoot circumference 18-19 cm)	#8185-18
Medium (fits midfoot circumference 20-21 cm)	#8185-20
Large (fits midfoot circumference 22-23 cm)	#8185-22

### Visco-GEL® Arch Support Wrap

Helps Relieve Heel and Arch Pain

Stretchable nylon fabric with hook and loop closures positions a soft silicone Gel cushion that supports the plantar fascia. Designed for wear under or over socks. Thin profile fits comfortably in most footwear. Washable and reusable for months. Interchangeable — fits right or left foot.



1/pk	Item#
Small/Medium (fits narrow and average width feet)	#1291-S
Large/XLarge (fits wide feet)	#1291-L



**Deluxe Metatarsal Cushion**

2 soft foam layers, hand-sewn nylon cover for added comfort, durability & reduced friction.

1/pk	Item#
Small	#8188-S
Large	#8188-L

interchangeable for right or left foot

**Forefoot Protection, Cushioning, Alignment, Support and Stability**

- Products for 'on-foot' and 'in-shoe' applications
- Relieve pain, improve comfort, and enhance foot function
- Absorb impact shock, shear and friction forces
- Provide additional plantar padding
- Ease metatarsalgia, fat pad atrophy, calluses, neuroma pain
- Reduce stress on metatarsal heads and shafts



**Visco-GEL® Metatarsal Pad with Toe Loop**

1/4" thick, All-Gel construction, stretchable toe loop.

1/pk	Item#
Small Left	#3005
Small Right	#3004
Large Left	#3007
Large Right	#3006



**Visco-GEL® Forefoot Cushion with Toe Strap**

1/8" Gel with top and bottom Soft'N Stretch fabric cover, elastic toe loop.


2/pk	Item#
Small	#1250
Large	#1252

interchangeable for right or left foot

**Economy Metatarsal Cushions**

Single layer foam with one-sided fabric cover for durability interchangeable for right or left foot

10/pk	Item#
OSFM	#8187



**Visco-GEL® Thin Forefoot Pad**

1/16" Gel with top and bottom Soft'N Stretch™ fabric cover.

2/pk	Item#
Small	#1275
Large	#1276

interchangeable for right or left foot



Toe Post Protection

#1131



#1133



Stops Chafing Between Toes

**Visco-GEL® Thong Sandal Guards**

Protect inter-digital space from sandal thong friction, pressure. Make sandals more comfortable. Self-adhering.

1 pr/pk	Item#
OSFM (Toe Protection Only)	#1131
OSFM (with Ball-of-Foot Cushion)	#1133



**Tuli's® Gel Metatarsal Cushions**

Patented 'waffle' design compresses, rebounds with each step.

1 pr/pk	Item#
OSFM	#P89



Contoured Pad Lifts Metatarsal Heads to Relieve Pain and Calluses

**FELTastic® Metatarsal Pads**

Soft support and cushioning. Self-adhesive backings. Ideal for neuroma pain, metatarsalgia, fat pad atrophy, 'dropped' met. heads, more.

1 pr/pk	Item#
Small (1/4" thick)	#8184-S
Medium (5/16" thick)	#8184-M
Large (3/8" thick)	#8184-L



Side View

**Self-Adhesive**



**New!**

**Pedi-GEL® Metatarsal Support Pads**

*Relieves Ball-of-Foot Pain, Calluses*

These soft, contoured Gel pads provide two benefits — they cushion and protect the forefoot from shock and pressure, and support the metatarsal or 'transverse' arch to help relieve calluses and neuroma pain. Self-adhering to stay in position in shoes. One size fits most.

1 pr/pk	Item#
OSFM	#8208





**New!**

**Reusable Gel Metatarsal Pads**

Soft, long-lasting Gel cushions adhere directly to the forefoot to relieve ball-of-foot pain, absorb, and redistribute pressure. Self-adhering. Washable, sturdy Gel lasts for months. Two thicknesses — specify Regular (1/8") or Thick (1/4"). Interchangeable, fits Right or Left foot.

2/pk	Item#
1/8"	#8218
1/4"	#8219



**New!**

**Reusable Gel Dancer's Pads**

Relieve forefoot and ball of foot pain by shifting weight away from the great toe joint and sesamoid bones. Adhere directly to feet. Use with or without shoes. Durable, washable Gel lasts for months. Natural stickiness is renewed with every washing. Two thicknesses – Regular (1/8") or Thick (1/4"). Specify Right or Left.

2pk	Item#
1/8" Left	#8214-L
1/8" Right	#8214-R
1/4" Left	#8215-L
1/4" Right	#8215-R



New Color Coming Soon!

**New!**

**Reusable Gel U-Shaped Callus Pads**

*Relieve, Remove, Wash, Reuse!*

Versatile horseshoe-shaped Gel cushions relieve calluses, ball-of-foot and forefoot pain, and sore spots on the feet and ankles. Naturally self-sticking material, adhesive renews with each washing. Effective for months.



2/pk	Item#
	#8213

**New!**

**Foam Metatarsal Pads**

Cushion, support, protect to stay in place and ease ball-of-foot and forefoot pain. Contains Latex. Self-adhesive. One Size.

8/pk	Item#
OSFM	#8603



**Ideal offloading components for PressureOFF® Insoles!**



**New!**

**Dancer's Pad,  
Premium Foam**

Soft, yet sturdy cushions shaped to off-load weight, pressure and impact from the big toe joint and sesamoids to relieve pain. Self-adhesive. Specify Right or Left.

4/pk	Item#
Left	#8604-L
Right	#8604-R



**New!**

**Dancer's Pads, Premium Felt**

Soft, yet durable FELTastic® pads ease metatarsal pain by off-loading pressure from the forefoot and great toe joint. Treat turf toe, sesamoiditis, hallux limitus, metatarsalgia, more. Adhesive-backed. Specify Right or Left.

4/pk	Item#
Left	#8532-L
Right	#8532-R



**New!**

**PolyPads™  
Metatarsal Pads**

Specialized foam pads ease forefoot and ball-of-foot pain by off-loading weight from the forefoot. Fit easily in most shoes. Self-adhesive. Small fits shoe sizes 4-7, Medium fits 8-10, Large fits 11+.

4/pk	Item#
Small	#8663-S
Medium	#8663-M
Large	#8663-L



**New!**

**PolyPads™ Metatarsal Bar**

The Metatarsal Bar is shaped differently than typical 'tear-drop' shaped pads, helping redistribute pressure off a larger area, and supporting the forefoot arch. Self-adhesive. Small fits shoes sizes 4-7, Medium fits 8-10, Large fits 11+.

1 pr/pk	Item#
Small	#8662-S
Medium	#8662-M
Large	#8662-L



**New!**

**Visco-GEL® ExtraArches™**


*Add Instant Support to Any Shoe Style*

Anatomically-shaped, fabric-covered Visco-Gel® longitudinal arch pads provide instant support, comfort and foot pain relief. Quick help for patients suffering symptoms of excessive pronation, heel, arch and ball of foot pain. Customized fit and support in any shoe style — without taking up forefoot or rearfoot space. Self-adhering strip keeps them in place. Sizing is approximate.

1 pr/pk	Item#
Small	#P1140-S
Medium	#P1140-M
Large	#P1140-L

**Sizing Chart**

Small fits shoe sizes 4-7  
Medium fits 8-10  
Large fits 10+



**New!**

**PolyPads™ Arch Cookies**

*Instant Support in Any Shoes*

Soft, contoured foam arch supports adhere in shoes to add instant support. They fit in shoes without the bulk of orthotics or insoles. Fast relief from arch pain, flat feet and fallen arches in children and adults alike. Uniquely durable foam material. Self-adhesive. Small fits shoe sizes 4-7, Medium fits 8-10, Large fits 11+.

1 pr/pk	Item#
Small	#8664-S
Medium	#8664-M
Large	#8664-L



**New!**

**PolyPads™ Heel Cushions**

*Absorb Shock to Relieve Pain*

Unique, durable material cushions, protects and absorbs shock. Adds padding to thin-soled shoes. Make standing and walking more comfortable. Great for pain relief from heel spur syndrome and plantar fasciitis. Adhesive-backed to stick in place. Sizing is approximate. 1/4" thick. Small fits shoe sizes 4-7, Medium fits 8-10, Large fits 11+.

1 pr/pk	Item#
Small	#8665-S
Medium	#8665-M
Large	#8665-L



**Visco-GEL® Plantar Fasciitis Arch Support Sleeve**

*Relieves Heel and Arch Pain*

Soft fabric sleeve keeps contoured gel pad in position under the arch for instant support and long-lasting comfort. Helps relieve arch, heel and plantar fascia pain. Interchangeable to fit left or right foot. Effective with or without socks and shoes.

1/pk	Item#
Small/Medium	#1450-S/M
Large/XLarge	#1450-L/XL



**FasciaFIX® Plantar Fascia Stretching Sock**

*Helps Relieve Heel & Plantar Fasciitis Pain Faster*

Soft, sock-like foot covering with adjustable strap that keeps the foot stretched while resting or sleeping. Helps prevent microscopic tearing of plantar fascia ligament, speeds healing.

1/pk	Item#
Small	#P6045-S
Medium	#P6045-M
Large	#P6045-L

**Sizing Chart**

Small fits Calf Circumference less than 7¾"  
 Medium fits Calf Circumference 7¾"– 15¾"  
 Large fits Calf Circumference 15¾"



**New!**

**Plantar Fasciitis Relief Strips**

*Immediate Relief for Heel & Arch Pain*

- Proven, Fast Relief from Heel Pain
- Minimize Stress on Fascia Ligament
- Trimmable, Self-Adhesive, Easy to Use
- Thin Design, Fit in All Shoes
- OTC-packaged for easier dispensing, retail referral or patient direct order.

*New! Extended Treatment 21-piece Pack!*

7/pk	Item#	21/pk	Item#
	#6060-7		#6060-21

## Treating Heel Pain Just Got Easier



### Peel-Away™ Adjustable Heel Lift

Simple, Multi-Layer Pad Solves Common Problems

- Reduces strain to ease Achilles tendonitis pain
- Compensates for leg length discrepancies to balance hips, back
- Elevates calcaneus to ease plantar fascia strain, heel pain
- 3/8" total thickness includes 2 removable 1/8" layers to adjust height
- Soft, durable construction



1/pk	Item#
Small (fits women's 4-7)	#6582-S
Medium (fits women's 8-10, men's 6-8)	#6582-M
Large (fits women's 11+, men's 9+)	#6582-L



### Biomechanical Heel Wedges

Limit Pronation Without Bulky Orthotics

- Firm rubber postings help control calcaneal eversion
- Reduce pronation-related stress on plantar fascia, ankles, knees, hips
- Use laterally to decrease adduction and ease pain from medial compartment knee syndrome and osteoarthritis.



1 pr/pk	Item#
Small (fits women's 5-7)	#6589-S
Medium (fits women's 8-10, men's 6-8)	#6589-M
Large (fits women's 11-13, men's 9-11)	#6589-L
XLarge (fits women's 14+, men's 12+)	#6589-XL



### Plastic Heel Cups

Help Ease Heel Pain

Easy customization! Heat in warm water, then mold to heel.

- Firm plastic construction disperses heel strike shock
- Gently squeeze fat pad under the calcaneus to utilize natural cushioning, shock absorption
- Helps relieve heel, shin, lower back, Achilles and ankle injuries

6 pr/pk	Item#
Small—Blue Star (fits women's up to 7)	#6575-S
Large—Red Star (fits women's 8+, Mens 6+)	#6575-L



Cushioning 'waffle' Design



Tuli's Original

### Tuli's® Heel Cups

Compressible 'Waffle' Shock Absorbers for Sore Heels

Patented design compresses, rebounds with each step

Tuli's Original 1 pr/pk	Item#
Regular (for body weight up to 175lbs.)	#P84-REG
Large (for body weight over 175lbs.)	#P84-LRG


\*quantity discounts:  
Order 1-3 dz @ \$6.15/pr  
Order 3 dz or more @ \$5.20/pr



Tuli's Gel

Tuli's Gel 1 pr/pk	Item#
Regular (for body weight up to 175lbs.)	#P85-REG
Large (for body weight over 175lbs.)	#P85-LRG

\*quantity discounts:  
Order 1-3 dz @ \$7.15/pr  
Order 3 dz or more @ \$6.20/pr



**New!**

**LowShock™ Gel Heel Straights™**

*Varus/Valgus Alignment, Now in Shock-Absorbing Pedi-GEL®!*

Typically used to help limit calcaneal eversion and pronation to improve foot function, relieve heel and arch pain. ‘Wedge’ design improves biomechanics. Soft, durable, self-sticking Pedi-Gel® absorbs shock. Discrete under shoe linings, without taking up space in footwear.

1 pr/pk	Item#
Small (fits shoe sizes 4–7)	<b>#P5180-S</b>
Medium (fits shoe sizes 8–10)	<b>#P5180-M</b>
Large (fits shoe sizes 10+)	<b>#P5180-L</b>

**NatraCure® Cold Therapy Socks**

*Targeted Cold Therapy for Relief of Foot Aches & Pains*

- Relieves pain from heel spurs, plantar fasciitis, toe problems, more
- Eases discomfort on bottom of foot, top of foot and Achilles’ Tendon
- Soothes tired, sore, aching feet, strains and sprains
- Place removable chilled or frozen Gel packs where needed
- Effective after injury or surgery
- Immediate relief of burning, tingling and neuropathic pain

1 pr/pk	Item#
Small (fits shoe sizes 6-9½)	<b>#705</b>
Large (fits 10+)	<b>#707</b>

(Sizing is approximate)



**EZ Mornings™ Plantar Fasciitis Nighttime Stretching Splint**

Stubborn plantar fascia condition? Morning heel pain a problem for your patients? Speed the healing process and let them wake up without discomfort. This unique splint maintains the foot in a dorsi-flexed position to stretch the plantar fascia while patients sleep. Helps minimize fibrous tears that occur when patients first walk after periods of rest. If you haven't tried a NightSplint, you haven't tried everything to relieve plantar fasciitis. Also effective for Achilles tendonitis. Soft, comfortable construction, adjustable hook-n-loop fastener, non-slip sole. Fits Right or Left foot.

1/pk	Item#
Small (fits up to Women's shoe sizes 9, Men's 8)	<b>#6040-S</b>
Large (fits Women's shoe sizes 9 ½+, Men's 8 ½+)	<b>#6040-L</b>

HPCS Code: L4396



**Proven Relief of Heel Pain**

**Exclusive fabric top-covers make these more comfortable!**



**GelStep® Heel Cups with Soft Spur Spot**


Proven most effective in clinical trials. Dual durometer cushioning.

1 pr/pk	Covered	Uncovered
Small	#5050-C-S	#5050-S
Medium	#5050-C-M	#5050-M
Large	#5050-C-L	#5050-L
XLarge	#5050-C-XL	#5050-XL
XXLarge	#5050-C-XXL	#5050-XXL

HCPCS Code: L3170




**Treatment Tip:**  
Some patients prefer wearing these under hosiery, against the skin, for better positioning.



**GelStep® Heel Cups with Soft Center Spot**

1 pr/pk	Covered	Uncovered
Small	#5025-C-S	#5025-S
Medium	#5025-C-M	#5025-M
Large	#5025-C-L	#5025-L
XLarge	#5025-C-XL	#5025-XL
XXLarge	#5025-C-XXL	#5025-XXL

HCPCS Code: L3170



3° Post to minimize pronation, flat feet or some valgus knee conditions

**GelStep® Posted Heel Pads with Soft Spur Spot**

Proven most effective in clinical trials. Dual durometer cushioning.

1 pr/pk	Covered	Uncovered
Small	#5075-C-S	#5075-S
Medium	#5075-C-M	#5075-M
Large	#5075-C-L	#5075-L

HCPCS Code: L3485

**GelStep® Heel Cups & Pads Sizing**

Small fits Women's Shoe Sizes 4-7  
 Medium fits Men's sizes 6-8 ½, Women's 7 ½ - 10  
 Large fits Men's sizes 9-11 ½, Women's 10 ½ - 12  
 Extra Large fits Men's sizes 12-14 ½, Women's 12+  
 XX Large fits Men's sizes 15+



**New!**

**Double-Action Gel Heel Cushions**

*Removable Disc to Ease Heel Spur Pain*

These soft Gel Cushions protect heels with every step. The exclusive circular wave pattern absorbs shock and disperses pressure to ease pain and make walking, running or standing more comfortable. If you have a heel spur, remove the center disc to offload pressure or pain from sore area.

- Use with or without removable disc
- Durable gel construction
- Gel 'grips' — won't slip out of place
- Fit in all types of shoes


1 pr/pk	Item#
Ladies'	#P5040-L
Men's	#P5040-M

**Absorb Shock for Heel Pain Relief!**



*Soft off-center, true dual durometer zone minimizes shock*

**Removable Discs**



**Treatment Tip:**  
Pads fit under, or on top of shoe insoles

*Blue Zone under heel*



**GelStep® Heel Pads with Soft Spur Spot**

1 pr/pk	Covered	Uncovered
Small	#5100-C-S	#5100-S
Medium	#5100-C-M	#5100-M
Large	#5100-C-L	#5100-L

HCPCS Code: L3485



**Medical-Grade Silicone Absorbs Shock Best**

GelStep® Heel Cups, Pads, Insoles and Orthotics are all made of genuine, medical-grade silicone which is proven superior as a long-term shock absorber to foams, minimizing forces from heels and feet up the kinetic chain. Silicone absorbs and disperses shock quickly, rebounding for the next heel strike.

**GelStep® Heel Pads with Soft Center Spot**

1 pr/pk	Covered	Uncovered
Small	#5105-C-S	#5105-S
Medium	#5105-C-M	#5105-M
Large	#5105-C-L	#5105-L

HCPCS Code: L3485

## GelStep® Silicone Insoles & Orthotics

Exclusive fabric top-covers make these more comfortable!

### Maximum Shock Absorption, Cushioned Foot Protection

Soft support for any patient, especially geriatrics. Choose specific features including arch and metatarsal supports, cupped heels and more.

- Absorb shock, cushion and protect
- Single and true dual durometer cushioning — softer **Blue Zones** for high impact areas
- Clinically-proven biomechanical forms
- Will never flatten, deform or degrade
- With Exclusive anti-microbial top covers that absorb odor, prevent bacterial growth, minimize friction and ease shoe donning

#### Relief for...

Arch pain	Metatarsalgia
Forefoot pain	Fat pad atrophy
Pronation	Morton's neuroma
'Flat feet'	Haglund's deformity
Plantar fasciitis	Arthritis
Heel Spurs	Knee, hip, back pain
Heel Pain	General foot discomfort



### GelStep® 3/4 Length Insoles

Comfort, support and protect to ease heel, arch and forefoot pain. Mild metatarsal pad. Soft heel and met. **Blue Zones**

1 pr/pk	Item#
Small	#7000-C-S
Medium	#7000-C-M
Large	#7000-C-L
XLarge	#7000-C-XL
XXLarge	#7000-C-XXL



### GelStep® Thin Dress Shoe Insole with Low Wide Met Rise

Dress shoe insole with low, wide metatarsal pad redistributes forefoot pressure

1 pr/pk	Item#
Small	#7029-C-S
Medium	#7029-C-M
Large	#7029-C-L

HCPCS Code: L3060



Protection and comfort in limited space. Fit in most fashionable shoes. Toe-free design offers more freedom for digits.

### GelStep® Insole & Orthotic Sizing

Small fits Men's Shoe Sizes 4-5 ½, Women's 5-7  
 Medium fits Men's sizes 6-7 ½, Women's 7 ½ - 9 ½  
 Large fits Men's sizes 8-9, Women's 10-12  
 XLarge fits Men's sizes 9 ½ - 11  
 XXLarge fits Men's sizes 11 ½ - 12 ½  
 XXXLarge fits Men's sizes 13+




**GelStep® Full Length Replacement Insoles**

Comfort, support and protect to ease heel, arch and forefoot pain. Soft heel and metatarsal head **Blue Zones**

1 pr/pk	Item#
Small	#7050-C-S
Medium	#7050-C-M
Large	#7050-C-L
Xlarge	#7050-C-XL
XXLarge	#7050-C-XXL
XXXLarge	#7050-C-XXXL

HCPCS Code: L3060

Softer Density 'Blue Zones'  
Absorb Shock



**GelStep® Insole & Orthotic Sizing**

- Small fits Men's Shoe Sizes 4-5 ½, Women's 5-7
- Medium fits Men's sizes 6-7 ½, Women's 7 ½ - 9 ½
- Large fits Men's sizes 8-9, Women's 10-12
- XLarge fits Men's sizes 9 ½ - 11
- XXLarge fits Men's sizes 11 ½ - 12 ½
- XXXLarge fits Men's sizes 13+



1stRaythotics, 3-Layer Toe Separators, '4 in 1' Super Soft Toe Cushions, 50/50 'Soft Formula', Action Orthotics, AjustaLoop, Align-A-Heel, Arch Binders, Arch Cradles, Bunion Care, Bunion Relievers, Callus Control Cream, Cracks & Calluses Foot Repair Cream, Clotrimazolil, Customizables, Daysplint, Deep-Healing Foot Cream, Dexterity, Diabetic Defense, Dual-Action Bunion Fix, Enhanced Silicon, Extra Arches, EZ Mornings, FasciaFIX, FELTastic, FootReacher, FootSplints, Forefoot Compression Sleeve, FungaFiles, FungaSoap, Gel Ultimates, Halluxthotics, GelStep, Hammer Toe Wrap, Heel Straights, Heel-So-Smooth, HV Hinge, Little Toe Sleeves, LowShock, MicroLoop, Microscopically-Clean, NightSplint, NoBend, PediFix, Pedi-GEL, PediPlast, Pedi-Quick, Pedi-Smart, Peel-Away, Podiatrists' Choice, PolyPads, Reattachable, SeamLess, Sesamoid Relief Sleeve, ShoeZap, SkinStrips, Soaking Crystals, Soft'N Stretch, SoftSplint, SoreSpot, Springthotics, Stay-Put Toe Separators, Stay-Put Toe Spacers, Step-N-Smooth, Stretch Socks, Tea Tree Ultimates, the No Sore Feet symbol, Toe Separators, Toe Trainers, ToeBuddy, ToeCoach, ToeCrutch, ToeWraps, Treat their shoes while you treat their feet, Tubular Foam Toe Bandages and Visco-GEL are trademarks of PediFix, Inc.

B2Splint is a trademark of JD & SW, LLC, Darco is a registered trademark of Darco, Inc., Just For Toenails is a registered trademark of Dr. Robert Spalding, LiquiCell is a registered trademark of Birchwood Labs, Lycra is a registered trademark of E.I. du Pont de Nemours and Company, NatraCure and SmartGel are registered trademarks and divisions of PolyGel, LLC, Plastazote and Poron are trademarks of Rogers Foam Corporation, ToeSnug is a registered trademark of Evenup, LLC, Tuli's and TuliGel are trademarks of International Comfort Products, Inc., Velcro is a registered trademark of Velcro Industries.



**New!**


**DiabeticDefense® Orthotic Insoles**

*Soft Support & Skin Protection!*

Dual layer – 3.0mm Pink Foam Top, 2.0mm EVA base -- combine to provide gentle medial support and effective friction reduction. They're a first line defense and help to our many aging patients with risk factors for skin breakdown, and in need of foot cushioning and arch support. These comply with Diabetic Shoe Program official requirements. Trim to fit precisely.



1 pr/pk	Item#
Small fits Women 7-8, Men 5-6	<b>#2340-S</b>
Medium fits Women 9-10, Men 7-8	<b>#2340-M</b>
Large fits Women 11-12, Men 9-10	<b>#2340-L</b>
XL fits Women 13-14, Men 11-12	<b>#2340-XL</b>
XXL fits Men 13-14	<b>#2340-XXL</b>



**Proven plantar fasciitis relief**

**TL-2100 Preforms® Orthotics**

*Maximum Control in Minimum Space*

Super-thin, technologically-advanced rigid carbon graphite material. Fit in all footwear. Heat moldable. Adult sizes.

**Item #4860**

Deep Heel Cups, No Metatarsal Raise

Women's Sizes 4-11  
Men's Sizes 6-13




*Super thin material*





**Custom Results, Preforms® Savings**




**Arch Cradles®  
Soft Support  
Preforms® Orthotics**

Soft Longitudinal and Metatarsal Support. Fit in most footwear.

1 pr/pk	Item#
Women's 5-6	#P182-XS
Women's 7-8, Men's 6	#P182-S
Women's 9-10, Men's 7-8	#P182-M
Women's 11+, Men's 9-10	#P182-L
Men's 11-12	#P182-XL
Men's 13-14	#P182-XXL



Deep Heel Cup




This product contains natural rubber latex which *may* cause allergic reactions.

**Cork & Latex Preforms® Orthotics**

Flexible Support, Deep Heel Cups  
Requires underside grinding to fit in most footwear.

**Item #1562** With Metatarsal Raise  
Women's Wide Width Sizes 4-11  
Men's Standard Width Sizes 6-13



**Leather Preforms® Orthotic Shells**

Create your own customizations — adhere cork, foam, felt, other materials.

**Item #1599**  
Women's Sizes 4-11  
Men's Sizes 6-13




**Order Preforms® by Exact Shoe Size**

*"...due to the low cost of your orthotics, I'm able to pass considerable savings along to my patients, and also give valid treatment to patients who don't desire surgery."*  
Dr. W., South Hadley, MA

### Action Orthotics™ 3/4 and Full-Length Supports

Relieve Heel, Arch & Foot Pain

These unique orthotics provide effective, economically-priced support for immediate comfort and foot function improvement. Innovative design helps improve biomechanics, especially during toe-off to reduce stress on the 1st MPJ. Ideal for bunion sufferers. ¾ Length style fits atop existing shoe insoles.



1 pr/pk	3/4-Length Item#	Full-Length Item#
Small	#P5610-S	#P5620-S
Medium	#P5610-M	#P5620-M
Large	#P5610-L	#P5620-L
XLarge	#P5610-XL	#P5620-XL



**New!**



#### Action Orthotics Sizing

Small fits Youth Shoe Sizes 1-4, Women's Shoe Sizes 5-7, Men's Shoe Sizes 4-6  
 Medium fits Women's Shoe Sizes 8-10, Men's Shoe Sizes 7-9  
 Large fits Women's Shoe Sizes 11-13, Men's Shoe Sizes 10-12  
 Extra Large fits Men's Shoe Sizes 13-15

Due to worldwide shortages of foam materials, colors may vary from those shown here and on packages.

*The comfort of a soft insole, with the control of a semi-rigid orthotic*

**Active Orthotics™ 3/4 and Full-Length Orthotics**

**Moderate Foot Function Control**

2mm thick thermoplastic provides comfortable support. Vinyl top cover. Fit in most footwear. Heat spot-moldable. Adult sizes.



*Provide support and comfort in most footwear styles.*

- Moderate foot function control
- Dual-layer foam with fabric top covers are durable
- Heat moldable shell for spot modification
- Deep heel cup for calcaneal control



**Item #5300**

Women's Size/Men's Size: W5, W6, W7, W8/M6, W9/M7, W10/M8, W11/M9, W12/M10, W13/M11, M12, M13

Please specify shoe sizes when ordering

**Item #5350**

Women's Size/Men's Size: W5, W6, W7, W8/M6, W9/M7, W10/M8, W11/M9, W12/M10, W13/M11, M12, M13

Please specify shoe sizes when ordering



*Full-length to replace shoe insoles*



**Plantar fasciitis relief!**

Due to worldwide shortages of foam materials, orthotic colors are subject to change.

# Visco-GEL® Ankle Protection Sleeves

**Gel padding inside!**



## Visco-GEL® Achilles Protection Sleeve

Relief from 'pump bump' (Haglund's deformity), sore, fatigued or injured Achilles tendon

- Soft, comfortable Gel pad is in direct contact with skin to protect the posterior of the heel from friction, pressure
- Slim profile design fits in most footwear, including boots, skates, etc.
- Use with our Peel-Away™ Adjustable Heel Lift (see page 36) to take pressure off the strained Achilles.

1/pk	Item#
Small/Medium	#1400
Large/XLarge	#1401

## Comfort, Relief and Protection for Ankle, Achilles, Dorsum and Malleolar Conditions

These devices feature our uniquely soft, compressive Soft'N Stretch™ fabric and our proprietary Visco-GEL® that conforms to, and protects specific areas.



## Visco-GEL® Malleolar Protection Sleeve

Relief and protection for the ankle bone

- Gel pad is in direct contact with skin to absorb pressure, friction, shock and shear
- Slim design fits in most footwear, including boots, skates, braces, etc.

1/pk	Item#
OSFM	#1405

*Slim design fits in most footwear*



## Visco-GEL® Dorsum Protection Sleeve

Relief and protection for the top of the foot.

- Absorbs pressure, friction, shock and shear
- Effective in shoes, boots, skates, etc.

1/pk	Item#
Small/Medium	#1408
Large/XLarge	#1409

## Achilles/Dorsum/Malleolar Sleeves Sizing Chart:

S/M fits Women's shoe sizes to 9, Men's to 7  
L/XL fits Women's shoe sizes 9+, Men's 7+

## Visco-GEL® Achilles/Dorsum Two-Way Protection Sleeve

Relief and protection for back of heel and top of foot.

- Gel pads are in direct contact with skin to absorb pressure, friction, shock and shear.
- Slim design fits in most footwear, also boots, skates, etc.

1/pk	Item#
Small/Medium	#1412
Large/XLarge	#1413



Now available in Black!  
Discreet compression support for fashion conscious patients.




#6021

**Pedi-Smart® Compression Ankle**

Helps Control Edema After Injury or Surgery

- Provides graduated compression from metatarsals to calf to assure increased venous activity
- Eases pain, reduces swelling, speeds healing
- Lightweight elastic stretches comfortably for optimal fit
- Use conservatively following ankle sprains, post-op after foot or ankle surgery
- Provides Grade III compression (30-40mm Hg)



#6020

Measure Here

New Design Improvements!  
Less Dorsal Bunching,  
Improved Heel Fit

Helpful for Plantar Fasciitis Treatment, Too

1/pk	Beige	Black
#2 - fits ankle circumference 20-22 cm (8"-8.75")	#6020-#2	#6021-#2
#3 - fits ankle circumference 22-24 cm (8.75"-9.5")	#6020-#3	#6021-#3
#4 - fits ankle circumference 24-26 cm (9.5"-10.25")	#6020-#4	#6021-#4
#5 - fits ankle circumference 26-28 cm (10.25"-11")	#6020-#5	#6021-#5
#6 - fits ankle circumference 28-30 cm (11"-11.75")	#6020-#6	#6021-#6

Material content: 100% uncoated spandex (synthetic hypoallergenic elasthane).

**New!**

**FasciaFIX® Plantar Fasciitis Relief Sleeve**

Eases Heel & Arch Pain

Soft, elasticized ankle with unique knitting pattern provides heel and arch compression to minimize swelling, support fascia ligament and relieve pain. Tactile 'nodules' help stimulate circulation. Worn over or under socks. Daytime or nighttime use to ease morning discomfort. Add this to your PF treatment protocol. Also helpful for Achilles tendon conditions. Provides 25-35 mmHG all over compression, with 30 mmHG at arch.



1/pk	Item#
Small	#P6023-S
Medium	#P6023-M
Large	#P6023-L
XLarge	#P6023-XL

**Sizing Chart**

Small: Arch Circumference 5"- 7.5", shoe sizes W4 - 6.5  
 Medium: Arch Circumference 7.5"- 10", shoe sizes M6 - 9.5 or W7 - 10.5  
 Large: Arch Circumference 10"- 12.5", shoe sizes M10 - 13 or W11+  
 Extra Large: Arch Circumference 12.5"- 15", shoe sizes M13.5+



Multi-Directional Compression

## Soothing Visco-GEL® Relief & Protection for Heels, Ankles, More



#1325



Before

After

Vented for breathability



#608

Moisture Zone

### Visco-GEL® Heel-So-Smooth® Heel Sleeves

Soothe & Soften Rough, Dry, Cracked Heels

- Lined with Visco-GEL® that releases vitamin E enriched mineral oil onto the skin surface
- Soften dry heels
- Use alone, or with heel creams for intensive treatment

1 pr/pk	Item#
Original	#1325
Vented	#608

Counter-Top Display	Item#
6 Pairs	#D09-HSS



Gel-Lined! Soothe & Soften

### Gel Ultimates® Moisturizing Booties

Softer feet in hours!

These Visco-GEL® lined slipper-socks moisturize and soften hard, dry, rough skin — fast. Get smooth feet after just a few short treatments. Exclusive Gel that fully surrounds feet — not just the soles — gradually releases jojoba, avocado, olive and other natural oils onto the skin's surface. Terrycloth comfort, non-skid soles. Safe to wear around the house while feet get super soft! Washable and reusable for months. One size fits most.

1 pr/pk	Item#
OSFM	#108



Fully-coated or half-coated with moisturizing Visco-GEL®

### Visco-GEL® Mesh Body Sleeves

Gel-lined Soft'N Stretch™ mesh fabric cushions & protects, soothes, compresses, relieves pain.

- 3" diameter, 10" length, *Trimable*
- Fits feet, ankles, calves, knees, wrists, elbows, arms
- Protects in skates & boots, under braces & splints

1/pk	Item#
Fully-coated	#1203
Half-coated	#1209

# Diabetic Solutions for Everyday, Active and Oversized Needs

**All our socks are Latex-free!**

## SeamLess™ Everyday Socks

Safer and more comfortable for people with sensitive feet or diabetes.

- Hand-sewn toe seams eliminate irritating ridges
- Cross-stretch design conforms to feet without binding or bunching

1/pk (white)	Item#
Small	#P795-S
Medium/Large	#P795-M/L
XLarge	#P795-XL

94% Cotton, 5% Nylon, 1% Elastic

- Knit padding cushions impact
- CoolMax® lining keeps feet dry and comfortable
- Easy on and off
- Cross-stretch design won't bunch
- Buried toe seam



## SeamLess™ OverSized Socks

Same features as the Everyday socks, plus extra-wide tops stretch comfortably over enlarged feet, ankles and calves, casts and bandages, too.

Resilient fibers allow them to return to their original shape after laundering.

1/pk (white)	Item#
Small	#P798-S
Medium/Large	#P798-M/L
XLarge	#P798-XL

96% Cotton, 4% Elastic

Stretches to 18 inches



## SeamLess™ Sock Sizing Chart

Small fits Women's Shoe Sizes 5-7 ½, Men's 4-6 ½  
 Medium/Large fits Women's Sizes 8-10 ½, Men's 7-9 ½  
 Extra Large fits Women's Sizes 11-13 ½, Men's 10-12 ½

## Stretch Socks™

Fit easily over enlarged feet, ankles, calves, casts and bandages.

- Resilient, non-bunching construction
- Micro-Cool™ treatment process moves moisture away from skin surface

1/pk (white)	Item#
Medium	#P760-M
Large	#P760-L
XLarge	#P760-XL

90% High-Bulk Acrylic, 8% Nylon, 2% Lycra Elastic

## Stretch Socks™ Sizing Chart

Medium fits Women's Shoe Sizes 7-10, Men's 6-9  
 Large fits Women's Shoe Sizes 10-13, Men's 9-12  
 Extra Large fits Men's Shoe Sizes 13+



# 'Zap' Shoes Clean to Fight Fungus, Infections, Odor & for Better Foot Health

15 Minute  
UV Shoe  
Sanitizer

# ShoeZap®



## Kills fungus & bacteria that cause...

- Athlete's Foot
- Toenail Fungus
- Diabetic Infections
- Foot & Shoe Odor

## The Problem — Our Shoes Are Contaminated

Warm, dark and damp, shoes are the perfect breeding ground for germs, including...

Fungus

Bacteria

Viruses

Yeasts

Molds

Allergens

*T. rubrum*

*T. mentagrophytes*




*P. aurgenosa*

*C. albicans*

*Staphylococcus aureus*

*E. coli*

More



If your patient's feet or nails are, or have been, affected by any of the above, it's likely their shoes are contaminated with these pathogens. Germ-filled shoes can be a source of foot and nail infections.

Fortunately, now there is help for people with these common foot conditions, and everyone wanting to keep their shoes cleaner, more hygienic and smelling fresh.

## The Solution — ShoeZap®

- Sterilizes microbes that can attack, infect, and re-infect feet
- Disinfects & Deodorizes — kills bacteria that cause odors
- Keeps shoes microscopically-clean™, more hygienic and safer to wear

## Proven technology — improved!

ShoeZap® is the first and only UV Sanitizer with all these unique advantages...

- Faster — works in just 15 minutes
- Simpler — ONE size Sanitizer fits ALL shoe sizes & styles
- Reversible — to fully disinfect the heel area, too
- Rugged — engineered for years of frequent use

## How ShoeZap® Works

ShoeZap® emits UVC germicidal light (253.7 nm wavelength) that sterilizes micro-organisms. This same technology has been used worldwide for decades to sanitize air in hospitals, kill germs in drinking water, disinfect food preparation surfaces and more.

## Science & Research

Numerous studies demonstrate the sanitizing power of UVC light. Others prove its effectiveness in footwear. See research, studies, expert opinions and more on our website here: <http://www.pedifix.com/t-szstudies.aspx>

## ShoeZap® — Instead of Toxic Shoe Sprays

Shoe sprays have limited, short-term effect on germs. Worse, the chemicals that might kill microbes can leave a toxic residue in footwear. Re-activated by perspiration, this creates an unhealthy 'foot bath' feet are then trapped in all day. ShoeZap® is 'green technology' — no smelly fumes, aerosol dispersants or lingering chemicals in shoes. And over time, ShoeZap® costs less, too.



Works in ALL types of footwear!



## Two Ways to Help Your Patients to Better Foot Health with ShoeZap®

### 1. "Treat Their Shoes While You Treat Their Feet"™

It works so fast, you can 'Zap' patient's footwear during their appointments. Doing so helps them understand its simplicity and benefits, and encourages purchase for home use. Complimentary ShoeZap treatments are highly appreciated by patients, and a foot health benefit.

### 2. Recommend, Dispense or Prescribe ShoeZap®

ShoeZap is available for in-office patient purchase — a convenience for them, with profits for your practice — or via patient direct order at 1-800-PEDIFIX, www.pedifix.com, and by mail, fax or email. It may also be available in stores near your practice for retail referral in the future (it's being introduced first to healthcare professionals).

## 15 Minutes and 4 Easy Steps to Microscopically-Clean™ Shoes

**1** Place Sanitizers in Shoes

**2** Place Shoes & Sanitizers in UV-Proof Safety Bags

**3** Power 'On'

**4** Device Turns Off Automatically After 15-Minute Treatment



### Patient Education Materials

(FREE with ShoeZap® purchase!)



ShoeZap Poster/Brochure Holder,  
20 Patient Brochures #LH2-SZ



ShoeZap Brochure Holder,  
with 20 Patient Brochures, #LH1-SZ



ShoeZap Patient  
Brochures,  
#LP3410-P

## Frequently Asked Questions

### Who should use ShoeZap®?

Everyone! Especially people susceptible to, wanting to prevent, or being treated for odor, fungus, skin, toenail, diabetic or post-op infections, those with compromised immunity, and anyone who wants to be proactive about their shoe cleanliness, foot health and hygiene.

### How often should I use ShoeZap®?

The more often, the better. Daily or every-other-day use is recommended to patients who wear the same shoes every day, don't wear socks, perspire heavily or have an active skin or nail infection. Weekly, monthly and seasonal use is suggested for shoes worn less frequently.

### How does ShoeZap® complement medical treatments for toenail fungus, Athlete's Foot and other infections?

Helping reduce the chance of recurrence and re-infection from contaminated shoes, ShoeZap® can enhance results from laser, oral and topical treatments.

ShoeZap® is the better, faster, safer way to kill germs that can infect your patient's footwear. And yours. Get ShoeZap for yourself, and your family, in addition to your patients and practice. Order today to reduce the chance of infection, and enjoy cleaner, more hygienic shoes.

To order, or get more information, contact us or your PediFix supplier today!



Description	Item#
ShoeZap® UltraViolet Shoe Sanitizer	#P3410
ShoeZap® UltraViolet Shoe Sanitizer, 10-Pack	#P3410/10
ShoeZap® Replacement UVC Bulb (Ozone-Free)	#P3411
ShoeZap® Replacement UV Protection Bag	#P3412

Each ShoeZap® Kit includes:  
2 UV Shoe Sanitizers (1 pair),  
2 UV-Proof Safety Bags,  
1 A/C Power Adapter, 1  
User Manual,  
1 Storage Box

*"Any treatment directed solely at the feet is inadequate... effective control measures must also include simultaneous eradication of the organisms in shoes"*

Journal of Investigative Dermatology

## Fungus? Here's Help! ClotrimazOil®

*The first and only broad-spectrum liquid antifungal that delivers Clotrimazole 1% in a proprietary, Essential Oils carrier formula*

*Frustrated with poor results from 'yesterday's' liquid antifungals — Tolnaftate & Undecylenic Acid — and the exorbitant cost of the latest 'wonder drug'? Here's a new alternative...*

Clotrimazole is active against 6 different causes of infection and discoloration, including fungus, mold and yeast. Tolnaftate and Undecylenic Acid topicals are not, perhaps explaining their limited effectiveness.

Those old-fashioned drugs are also frequently delivered in water or alcohol carriers, neither providing any inherent or additional antifungal benefits. Our scientists helped pioneer a new process to make Clotrimazole soluble in Essential Oils with a history of antimicrobial, antiseptic, antifungal, conditioning and penetrating activity. The resulting combination is ClotrimazOil®. It is now recommended widely as a primary antifungal topical therapy, a complement to oral drugs and laser nail treatments, and a long-term recurrence preventative. It is FDA-approved, safe and effective.

### No Reason to be 'Jubilant' — Still No Cure for Fully-Involved Onychomycosis

"If TV ads and free steak dinners correlated with improved treatment outcomes, we'd all be fat and rich" said my colleague. But no matter how many meals 'big pharma' buys us, or TV ads they run, there is still no cure for fungal nails we see in our practices every day.

It's hard to get excited about the latest 'wonder drug' — with complete cure rates under 18%, after 48 weeks of use, on 'mild to moderate' infected nails only 20-50% involved. How many of your patients seek treatment at that low infection threshold? Probably very few.

At a cost of hundreds of dollars per tiny tube and thousands more over the one year treatment time — generally NOT paid for by insurance — there is little to celebrate in those results, especially considering the likelihood of recurrence and re-infection.

Instead, improving appearance, symptom management and reducing the fungal burden are realistic outcomes for most patients.

ClotrimazOil® is a more cost-effective solution, and the first and only broad-spectrum 'azole' antifungal in this proprietary carrier formula.

So go ahead, enjoy a free steak dinner. But read the fine print in the study results. There's no reason to be jubilant. Switch to ClotrimazOil® — from old-fashioned Tolnaftate and Undecylenic Acid-based topical antifungals or over-priced 'wonder drugs', neither proving effective for our typical patients.

### Exclusive Carrier Formula

ClotrimazOil® delivers Clotrimazole 1% to affected surfaces in a proprietary Essential Oils blend that includes Tea Tree, Jojoba, Peppermint, Sweet Almond and others.

#### Tea Tree Oil

Effective in vitro against dermatophytes found on the skin. In studies, it has been shown to stop bacteria, mold and fungal growth.



#### Jojoba Oil

Known as a fungicide, it is very similar to the sebum oil produced naturally by human glands, resulting in its popularity in many skin-care products.



#### Peppermint Oil

Has a long tradition of medicinal use, and is thought to activate cold-sensitive receptors in skin tissue, which is the primary source of the cooling sensation that follows its topical application.



### Frequently Asked Questions

#### How should my patients use ClotrimazOil®?

By applying a thin coating to affected surfaces — skin between toes, on, around and under nails, wherever reachable with the applicator brush — twice daily — morning and night. Use sparingly — this is a concentrated formula. Allow to dry.

#### How long will it take to work?

Fungal infections are difficult to eliminate. Patients should notice results within days on infections between toes, with continued improvement over time. If used on, around and under nails, it generally takes about 1 year for nails to re-grow. Encourage your patients to keep using ClotrimazOil® even after their symptoms are gone — fungus will return if not thoroughly treated.

#### Can children and pregnant women use ClotrimazOil®?

Clotrimazole is safety rated 'B' for both groups, unlike other topical antifungals.

#### How long will a bottle of ClotrimazOil® last?

3-6 months or more, depending on the amount used daily and the size of the affected areas it is applied to.

#### What if my patients aren't satisfied?

PediFix offers a Patient-Direct Satisfaction Guarantee. If for any reason you or your patients are not satisfied with the performance of ClotrimazOil®, we will refund the cost of their purchase, directly to them or you. For more information, call 1-800-424-5561.



*"The most effective topical antifungal formulation I've ever used, and well worth a few extra pennies per day"*

Dr. C., Deerfield Beach, FL

*"I've used Tolnaftate products for so many years and they rarely work. I'm expecting superior clinical results with ClotrimazOil"*

Dr. H, Westchester, NY

*"Double-blind, multicenter, randomized, controlled trial of 117 patients found cure rates of 11 percent in the clotrimazole treated group and 18 percent in the tea tree oil group. Clinical assessment showed partial to full resolution in approximately 60 percent of the patients in both treatment groups."*

Buck DS, Nidorf DM, Addino JG. Comparison of two topical preparations for the treatment of onychomycosis: Melaleuca alternifolia (tea tree) oil and clotrimazole. J Fam Pract. 1994;38(6):601-605.



**New!**

**ClotrimazOil®  
Topical Antifungal Solution**

*The first topical with the broad spectrum antifungal power of Clotrimazole, in an exclusive penetrating carrier oils formula.*

Highly effective on infected surfaces between toes, around, adjacent to and under nails where reachable with the convenient applicator brush.

- Kills 6 types of fungus, molds and yeasts
- Exclusive Essential and Carrier Oils formula
- Natural antioxidants and skin conditioners
- Urea to help soften nails
- Paraben-Free
- FDA-approved lab
- Made in the USA
- Doctor Formulated
- Patient-Direct, Money-Back Satisfaction Guarantee

.5 Oz. Bottle	Item#
1 bottle/box	#3490
12 bottles/case pack	#3490/12

**New Penetration Technology — Terpene Rich Essential Oils**

**Exclusive Formula with Tea Tree Oil!**

**Effective on webspace infections, too**

**Use with FungaFiles® to help remove fungal nail debris see page 61**





**New!**  
**Just for Toenails®**

**Enhanced Nail Polish**  
**Tea Tree Oil Enriched**

- Improve the Appearance of Problem Nails
- 10-Chemical Free
- Vegan Friendly
- O<sub>2</sub> Permeable
- Podiatrist Developed
- Healthy Nail Formula



New Bottle Design Coming Soon!

Loose Lips Pink Ships — #484	I Love You From My Head Toe-Ma-Toes — #569	Fuschia Flash — #121	Robin The Cradle — #430	Shaken Not Stirred — #571	Teal Of The Century — #479
Smoky Rose Mauve — #339	Private School Girl — #359	Plant Your Sugar Lips Here — #520	Pelican Peak Purple — #522	Paint Me Now and Later — #523	Icing On The Cake — #560
Burnt Orange — #600	Royal Fool's Gold — #313	Hickory Trickery — #420	Concrete Jungle — #354	Smoke Signal — #410	Good For Muffin — #400
Fiery Coral — #324	Lucky in Atlantic City — #491	Mad Money — #380	Mountain View Blue — #384	Orange You Sexy — #284	Clear Coat — #200
South Beach Pink — #357	Sub-Lime — #360	Sunshine & Daisies — #442	Tiramisu — #254	White — #280	Base Coat — #100

**Durable Colors for Pretty Toes & Feet All Year**

A renowned Podiatrist just couldn't take it any longer! Patients arrived in his office with unsightly, fungal nails, signs of infection from poor salon hygiene, and chemical-filled polishes that only made nail problems worse. After years of frustration, he developed his own exclusive, healthy nail polishes, with nearly instant results — improvements in nail appearance, reduced discoloration, normal thickness, less brittle, too! Discover a solution to a common, embarrassing problem. Free of Formaldehyde, Phthalate, Toluene.

1/2 ounce bottles. Specify colors by #.

**1 bottle Item#**  
 0.5 oz. **#M905**

**Ask about display options for your practice!**

Colors Change Frequently!

[www.pedifix.com/t-JustForToenails.aspx](http://www.pedifix.com/t-JustForToenails.aspx)

♥ Most Popular



**5 1/2"**  
**Professional Nail Cutter**

*Trims Hard, Thick, and Fungal-Infected Nails with Ease*

- 'Double-Action' Design For Extra Nail Cutting Power
- Heavy-Duty, Stainless Steel
- Etched Sure Grip Handle
- Same Design Podiatrists Use
- Packaged for Sale to Appropriate Patients for At-Home Nail Care
- Tough Enough for Clinical Use

1/pk	Item#
Curved	#P3355-CRV
Straight	#P3355-STR



**Choose Curved or Straight Blade Styles!**




**38% Urea Formula!**

**New!**  
**Urea Nail Butter**

*Urea Stick Dissolves Fungal Nails*

Lavender oil, tee tree oil and Urea combine to improve the appearance of thick, yellow toenails, safely and effectively. .19 oz.

2/pk	Item#
	#3062



**FungaFiles®**  
*First Disposable Files for Fungal Nails & Warts*

Encourage patient at-home debridement between visits to speed treatment results. Sized for precise use on varying surfaces these extra-gritty, double-sided, miniature files help remove fungal debris to help nails look better and make them thinner for easier trimming. Effective on warts to remove dead tissue, allow topical agents to penetrate deeper. Unique 24-piece, 'matchbook' package encourages frequent use and disposal. Complements oral, topical and laser therapies.



12 packs/bag	Item#
	#P3012-12

**Use with ClotrimazOil®**  
see page 59

# Relieve & Prevent Stubborn Skin Conditions

## Tea Tree Ultimates® FungaSoap®

*Cleans away fungus and bacteria*

- Studies show that Tea Tree Oil kills fungus and bacteria that can contribute to skin infections\*
- Regular washing with FungaSoap® may significantly reduce the occurrence and recurrence of infections
- Results from thousands of users suggest high efficacy, high compliance and no side effects



Convenient for use in bath or shower

\*Studies available on request... <http://www.pedifix.com/t-FungaSoapStudies.aspx>

1 Bottle	Item#
6 oz.	#P3071
13.5 oz.	#P3072
2 oz.	#P3075

*"It really works. Following treatment, I tell my patients to plan on using this for life."*

— Dr. José Concha

**New!**

## FungaSoap® CBD Relief Formula

The addition of 200 mg of U.S. grown and processed Full Spectrum CBD enhances the FungaSoap experience! Combined with the effects of Tea Tree Oil, this invigorating foot wash leaves feet feeling refreshed, clean, and smelling great.

1 Bottle	Item#
6 oz	#P3085



## Diabetic Defense® Daily Therapy Foot Wash

*Extra Protection for At-Risk Feet*

Physician-formulated with 12 specialized and Essential Oils. Tea Tree Oil cleans and protects. Ginger stimulates circulation. Shea Butter with vitamins A, E & F to moisturize. Shower, bath or tub convenience. Regular use helps soothe, soften, moisturize skin and helps minimize infection risk by washing fungus and bacteria from the skin surface.

*With Ginger Oil to help improve circulation, and Shea Butter for extra moisturizing*

1 Bottle	Item#
5.1 oz.	#P3076
Counter-Top Display	Item#
9 – 5.1 oz. bottles	#D09-DDFW

**New!**

## Dermatology 'Trial Size' Sample Packets



100 Packets	Item#
5mL Deep-Healing Foot Cream™	#P3069-PACKET
3mL FungaSoap®	#P3071-PACKET
5mL Cracks & Calluses Foot Repair Cream™	#P3305-PACKET
5mL Callus Control Cream™	#P3310-PACKET

**Exclusive Salicylic & Lactic Formula!****New!****Cracks & Calluses Foot Repair Cream™***Relieves Stubborn Skin Conditions*

Custom-compounded formula with Salicylic and Lactic acids in a petrolatum base, proven extremely effective in long-term clinical use for relieving rough, dry, callused, cracked and fissured feet and heels.

**1 Jar**    **Item#**  
 4 oz.    **#P3305**

**Exclusive Vitamin Formula!****Deep-Healing Foot Cream™***Softens rough, dry, cracked skin*

- Moisturizes calluses and heel fissures
- Penetrates fast and deep
- Non-greasy formula with vitamins A, D and E

**1 Jar**    **Item#**  
 4 oz.    **#P3069**

**Exclusive Shea Butter Formula!****Diabetic Defense® Daily Therapy Foot Moisturizer***Developed for Patients with Diabetes, But Everyone Loves It!*

Exclusive formula penetrates deep to soften thick calluses and moisturize dry, cracked skin and fissures. Contains Shea Butter with vitamins A, E, & F for long lasting hydration. Non-greasy, absorbs quickly.

**1 Jar**    **Item#**  
 4 oz.    **#P3620**

**Exclusive 20% Urea Formula!****Podiatrists' Choice® Callus Control Cream™***20% Urea to Soften & Moisturize Rough, Dry Skin*

Hydrate and exfoliate in one step with this extra-strength formula. Urea penetrates calluses and thick, hard skin to produce soft, smooth feet that look and feel healthy.

**1 Jar**    **Item#**  
 4 oz.    **#P3310**



**Treatment Tip:**  
Use ½ packet for foot soaks



**Tea Tree Ultimates® Soaking Crystals™**

*Therapeutic foot soak that softens, soothes, helps protect skin from fungus & bacteria*

- Exclusive formulation of Epsom and Sea Salts, Tea Tree and Peppermint essential oils
- Helps cleanse and soothe to maintain soft, healthy skin
- Reduces inflammation, fights infection, relieves odor, itching and pain
- Helpful for multiple skin conditions and after nail avulsion to promote healing and reduce the chance of recurrence

1 Box	Item#
(6) 1 oz pkts	<b>#P3078</b>
Bulk Pack	Item#
(100) 1 oz pkts	<b>#P3078/100</b>



Counter-Top Display	Item#
6 Boxes	<b>#D09-SC</b>

#D09-SC



**Tea Tree Ultimates® Foot Wipes**

*Deodorize & Refresh Feet, Fast!*

- Wipe Away Sweat & Odor, Dirt & Grime, Fungus, Bacteria & Germs
- Exclusive Footcare Formula: Tea Tree & Peppermint Essential Oils, Aloe & Vitamin E — Alcohol Free
- Pre-exam foot cleaning by doctor, staff or patient — keep a packet in every treatment room
- Quick, easy foot cleanings in hospitals, nursing homes, long-term care facilities
- Diabetic patient foot care & education — ‘Clean, Inspect, Prevent’
- Post-op or wound care around incision sites, ulcers
- Whenever a shower is unavailable or inconvenient

25 Wipes/pk	Item#
6.6" x 7.8" (16x20cm)	<b>#P3370</b>

24 packs/case	Item#
	<b>#P3370/24</b>



*“Keeping feet clean is key to foot health. These new Wipes are a big help to my patients.”*


R. Maldonado, DPM



**SoreSpot® Advanced Blister & Skin Protection Bandages with LiquiCell® Technology**

Unique liquid-filled Bandages absorb friction and shear forces. Cushion and protect to help prevent or relieve blisters and skin injury. Flexible fluid chamber provides a dynamic, shape-changing, bursa sac-like shield. Soothes corns, calluses, bunions, heel sores, diabetic 'hot spots', more. Usually applied to skin, but also effective on insoles/orthotics, even inside shoes. Clear, virtually invisible. Self-adhesive.

Size	Item#
Small, 4/pk	#P810-S
Large, 4/pk	#P810-L
Small, 10/pk	#810-S-10
Large, 10/pk	#810-L-10
Small, 100/pk	#810-S-100
Large, 100/pk	#810-L-100



**New!**

**Semmes-Weinstein Sensory Assessment Tool**

Map Sensory Changes, Detect Potential Trouble Spots Early


Test each patient's sensation response — quickly, easily — to diagnose nerve or neuropathy conditions in advance of complications. 10 gram monofilament. Individually packed with Testing Site Guide.

20/Box Item# #9005-20



Super Thin

Liquid-Filled Chamber 'Floats' Friction from Skin Surface



**Pedi-Quick® Dual-Action Foot File**

Quickly Reduces Corns and Calluses

Stainless steel 'coarse' surface on one side for faster, more aggressive action, 'fine' emery surface on the other, for gentle filing.

1/pk Item# #P3011



**Pedi-Quick® 2-Sided Foot File**

Works Gently, Safely

Two different emery surfaces — 'fine' and 'coarse' to smooth rough, dry, hardened skin, corn and calluses.

1/pk Item# #P3020

**German Made**



**Pedi-Quick® Step-N-Smooth® Shower Foot File**

Smoother Feet — No Bending!

Simply rub feet against this unique surface to safely and easily smooth calluses and rough, dry skin. Long-lasting mat rests on shower floor. Non-skid backing keeps it securely positioned. Great for people who cannot easily bend to reach their feet.

1/pk Item# #P3030



## OTC-Packaged PediFix® Footcare Products ... Recommend Them with Confidence

The products you trust, OTC-packaged, for easier, more profitable dispensing from your office or clinic, referral to a nearby store, or patient direct ordering. You choose how best to get PediFix products to your patients.



4' Medical Footcare Center  
#D40

*"I know your quality is unmatched, so PediFix is the only brand I suggest to my patients."*

Dr. C. Phoenix, AZ

## Dispense

### Dispense PediFix® OTC-Packaged Footcare Products from your office or clinic...

- Patients appreciate the convenience of getting these from you,
- Attractively packaged for easier sale,
- Each package offers directions for use, care and warnings to increase patient compliance,
- They're easier to charge for, reducing your costs of giving away treatment products for free (instead, generating sales and profits for your practice),
- Consider Dispensary Displays organize, stock and show these items in your office or clinic.



4-Sided Counter-Top  
Spinner Display  
#D27

## Prescribe

### Prescribe...through a nearby, locally owned retailer, not big box chains or Amazon — too many counterfeit products that will disappoint

- Get Prescription Pads, Patient Product Brochures, samples, more.
- <https://stores.pedifix.com/StoreFinder>, or call us to locate or arrange local availability in an independent pharmacy, DME or comfort shoe store near you.



Prescription Pads  
#L2041

## Patient Direct Order

### Enable your patients to purchase PediFix® OTC-Packaged products directly from PediFix...

- Call for Direct Order Brochures you can give to patients
- Guide patients to the consumer/patient portal at [www.pedifix.com](http://www.pedifix.com).



Product Guide  
#LPG—  
Available with or  
without Patient  
Direct Order Form

[www.pedifix.com](http://www.pedifix.com)  
New & Improved!

**2-Sided Spinner  
Counter-Top Display**



#D17

This two-sided spinner display fits easily on a counter, table-top or deep shelf to offer your patients relief and comfort from common foot ailments with the most popular PediFix medical footcare products.

**#D17** —  
2-sided counter-top spinner display  
22" H x 13.5" W x 12" D

- Display Includes:
- 20 Products (3 pieces each, assorted sizes as applicable)
  - 1 double-sided Display Header (16" x 9")
  - 25 Product Guides/Referral Brochures & Holder
  - Ailment Guide Shelf Hanger with lanyard
  - Out of Stock Tags

Product assortment and price subject to change. Request a current contents list today. Customize a display assortment with your favorite products!  
See Contents List at:  
[www.pedifix.com/t-displayinfo.aspx](http://www.pedifix.com/t-displayinfo.aspx)

**In-Office Dispensing  
Simplified & Professionalized  
with PediFix® OTC**

Giving treatment products away for free? Sending patients online to find them? Instead, boost compliance by offering the medical footcare products they need, at the time of their visit. These professional displays transform small space — in waiting or treatment rooms, on or near a reception desk or nurse's station — into a productive footcare product center.

- Convenience for your patients...boosts compliance, for better outcomes
- Revenue for your practice...easier to charge for
- Professional packaging enhances your image

**Choose a display for your office or clinic today!**

**Floor Spinner Display**



#D77

*"I rely on PediFix products for many of my patients. In-office dispensing makes it so easy to be sure that they get the right products, right away."*

Dr. G., Tamarac, FL

This mobile, two-sided 'spinner' features 50 of the most popular PediFix® products, in a space-saving 'footprint'. Includes a small assortment of insoles, heel cups and skincare items on bottom shelves. Heavy-duty, quality construction.

**#D77** —  
Floor Spinner Display  
72" H (w/header), 16" x 19" base with wheels

- Display Includes:
- 54 Products (4-6 pieces each, assorted sizes as applicable)
  - 1 double-sided Display Header (16" x 9")
  - 1 Window Decal
  - 25 Product Guides/Referral Brochures & Holder
  - Ailment Guide Shelf Hanger with lanyard
  - Out of Stock Tags

**Guaranteed success  
in your practice or  
we buy it back!**



OTC #P810, S, L, 4/pk  
Professional #810, Page 65



OTC #P81, S, L, 2/pk  
Professional #1000, Page 25



OTC #P30, M, L, 2/pk  
Professional #1006, Page 26



OTC #P32, 2/pk  
Professional #1026, Page 25



OTC #P338-MIX, 2/pk  
Professional #1027, Page 25



OTC #P27, M, L, 2/pk  
Professional #1031, Page 27



OTC #P33, 1/pk  
Professional #1033, Page 26



OTC #P35, 1/pk  
Professional #1035, Page 26



OTC #P82, S, L, XL, 1/pk  
Professional #1050, Page 23



OTC #P80, S, L, XL, 1/pk  
Professional #1075, Page 22



OTC #P31, S, L, XL, 2/pk  
Professional #1100, Page 24




OTC #P28, S,M,L, 2/pk  
Professional #1125, Page 27




OTC #P1133, pair  
Professional #1133, Page 37



OTC #P1202, S, L, 1/pk  
Professional #1202, Page 25




OTC #P1203, 1/pk  
Professional #1203, Page 54




OTC #P1225, SL, SR, LL, LR, 1/pk  
Professional #1225, Page 34



OTC #P1291, S, L, 1/pk  
Professional #1291, Page 25



OTC #P1295, S, L, 1/pk  
Professional #1295, Page 34



OTC #P1303, S, L, 1/pk  
Professional #1303, Page 31




OTC #P47, 1/pk  
Professional #1315, Page 30



OTC #P46-L, 1/pk  
Professional #1316, Page 30



OTC #P46-S, 1/pk  
Professional #1319, Page 30



OTC #P1321, SL, SR, LL, LR, 1/pk  
Professional #1321, Page 34



OTC #P1400, S, L, 1/pk  
Professional #1400, Page 52





OTC #P1405, 1/pk  
Professional #1405, Page 52



OTC #P1408 S, L, 1/pk  
Professional #1408, Page 52




OTC #P3095, roll  
Professional #3095, Page 12



OTC #P3004, SL, SR, LL, LR, 1/pk  
Professional #3004, Page 36



OTC #P53, S, M, L, Left or Right  
Universal, OSFM, 1/pk  
Professional #3025, Page 20



OTC #P805, /pk  
Professional #3075, Page 14



OTC #P806, 2/pk  
Professional #3076, Page 14



OTC #P60-OSFM, 1/pk  
OTC #P60-2, S, M, L, XL, 2/pk  
Professional #6000, Page 35



OTC #P6020, S, M, L, XL, XXL, 1/pk  
Professional #6020, Page 53



OTC #P6035, S, M, L, Right  
or Left, 1/pk  
Professional #6035, Page 32



OTC #P6040, S, L, 1/pk  
Professional #6040, Page 43



OTC #P51, 2/pk  
Professional #6051, Page 21




OTC #P6582, S, M, L, 1/pk  
Professional #6582, Page 42



OTC #P316, S, M, L, pair  
Professional #6589, Page 42



OTC #P15, 10/pk  
Professional #8005, Page 17



OTC #P16, 8/pk  
Professional #8023, Page 17



OTC #P229, S, M, L, 2/pk  
Professional #8129, Page 27



OTC #P280, S, M, L, Mix, 6/pk  
Professional #8130, Page 24



OTC #P34, S, M, L, 1/pk  
Professional #8134, Page 22



OTC #P337, S, M, L, Mix, 3/pk  
Professional #8137, Page 24



OTC #P54, S, M, L, XL, Left or Right, 1/pk  
Professional #8154B, Page 20



OTC #P55, 1/pk  
Professional #8155, Page 21




OTC #P57, 1/pk  
Professional #8157, Page 21



OTC #P58, Right or Left, 1/pk  
Professional #8158, Page 21



OTC #P86, S, M, pair  
Professional #8184, Page 37



OTC #P88, 1/pk  
Professional #8188, Page 36



OTC #P83, 2/pk  
Professional #8230, Page 24



OTC #P8201, 2/pk  
Professional #8201, Page 15




OTC #P8202, 2/pk  
Professional #8202, Page 15




OTC #P8203, 2/pk  
Professional #8203, Page 15




OTC #P8204, 2/pk  
Professional #8204, Page 15



OTC #P8205, 2/pk  
Professional #8205, Page 14



OTC #P8206, 12/pk  
Professional #8206, Page 15



OTC #P61, S,L, 1/pk



OTC #P68, 2/pk



OTC #P90, pair





OTC #P78, S, M, L, 1/pk  
Professional #1078, Page 22




OTC #P2105 6/pk  
Professional #2105, Page 16



OTC #P8501, 2/pk  
Professional #8501, Page 16



OTC #P8502, 2/pk  
Professional #8502, Page 16



OTC #P8504, 8/pk  
Professional #8504, Page 16



OTC #P8530, S, M, L, 2/pk  
Professional #8530, Page 16



OTC #P8534, 2/pk  
Professional #8534, Page 16



OTC #P8536, 2/pk  
Professional #8536, Page 16



OTC #TBD, 1/pk  
Professional #TBD, Page 2



OTC #P4862, S, M, L, pair  
Professional #4862, Page 9





OTC #P231, S, M, L, XL, XXL, pair  
Professional #2301, Page 11



OTC #P230, S/M, L/XL, XXL, pair  
Professional #2300, Page 11



OTC #P5620, S,M,L,XXL pair  
Professional #P5620, Page 50



# Complete Medical Footcare Center



## Help Your Patients Help Themselves With Our Growing OTC Assortment

**Dispense** For the highest compliance and best outcomes, stock and dispense PediFix products in your office. Choose OTC or traditional medical packaging options.

**Prescribe** Refer patients to locally-owned retailers (not chain stores) to get the products they need, immediately. See our online store locator, or call for the names and locations of stores near you that carry PediFix.

**Patient Direct Order** Enable your patients to purchase directly from PediFix via phone 1-800-PEDIFIX, free Patient Brochures, or online at [pedifix.com](http://pedifix.com) (not amazon — too many counterfeit products).

No matter what your practice setting, give your patients the benefits of long-lasting relief. Help them 'get the good stuff' — PediFix Medical Footcare products.

### New!

**#D40-3** —  
Full Line Floor Display  
60" H x 144" W x 17" D

Includes:

- 146 Products (2–12 pieces each, assorted sizes as applicable)
- 3 double-sided Display Header (16" x 9")
- 1 Window Decal
- 25 Product Guides/Referral Brochures & Holder
- Ailment Guide Shelf Hanger with lanyard
- Out of Stock Tags

This, and smaller OTC displays, are practical additions to a foot & ankle clinic. See pages 65-67 or call for details.

\*Display contents are subject to change. We introduce new products every year so the merchandise assortments and exact display merchandise price may vary slightly from those included here. See Contents List at [www.pedifix.com/t-displayinfo.aspx](http://www.pedifix.com/t-displayinfo.aspx)

# Helpful New Designs



**New!**  
**Visco-GEL® 3-Loop Hammer Toe Crest Cushion**  
*Lifts Toe Tips, Engages Toes*  
 #1036

PAGE 20



**New!**  
**Halluxthotics™ Turf Toe Insoles**  
*Limit Big Toe Motion to Ease Pain*  
 #4862

PAGE 9



**New!**  
**1/4" Lateral Sole Wedge Insoles**  
*Beveled to Limit Supination*  
 #2300

PAGE 11



**New!**  
**FungaSoap CBD**  
*Tea Tree Oil and CBD Combination!*  
 #P3085

PAGE 62



**New!**  
**Visco-GEL® Sport Toe Cap**  
*Extra Thick for More Protection, Durability*  
 #1078

PAGE 22



**New!**  
**1stRaythotics™ Preforms®**  
*Integrated Reverse Morton's Extension and Arch Support Relieve Pain*  
 #4875

PAGE 6



**New!**  
**FELTastic® Arch Cushions**  
*Instant Support Makes Feet Feel Better*  
 #8536

PAGE 16



**New!**  
**FELTastic® Heel Pad**  
*Lifts & Cushions to Relieve Heel, Calf, Body Pain*  
 #8534

PAGE 16



**New!**  
**FELTastic® Oval Callus Pads**  
*Cushion & Protect Sore Spots*  
 #8504

PAGE 16



**New!**  
**1stRaythotics™ Preforms®**  
*Metatarsal/Neuroma Offloader Relieves Pressure & Pain*  
 #4874

PAGE 7



**New!**  
**FELTastic® Metatarsal Pad, Contoured**  
*Supports Forefoot Arch*  
 #8530

PAGE 16



**New!**  
**Action Orthotics™**  
*Full-Length Affordable Support*  
 #P5620

PAGE 50



**New!**  
**Full-Length Active Orthotics**  
*High Performance Everyday Support*  
 #5350

PAGE 51



**New!**  
**Active Orthotics™**  
*3/4-Length Support, Control, Comfort*  
 #5300

PAGE 51



**New!**  
**1stRaythotics™ Preforms®**  
*High Performance Medial Support*  
 #4873

PAGE 7

**Clever New Ideas**



**New!**

**Bunion Comfort Splint™**  
Relieves, Aligns, Conservatively or Post-op  
#6025

PAGE 32



**New!**

**HVHinge Splint™**  
Aligns Big Toe, Allows Plantar  
& Dorsi Flexion  
#6038

PAGE 32



**New!**

**Moleskin-Lined Toe Separators**  
Divide, Align, Absorb Moisture  
#8131

PAGE 24



**New!**

**Urea Nail Butter**  
38% Urea Formula Dissolves Nails  
#3062

PAGE 61



**New!**

**1stRayometer™**  
Measure Range of Motion  
Accurately, Consistently  
#9000

PAGE 7



**New!**

**Neuropathy Monofilament Tester**  
Fast, Easy Assessment of Sensory Loss  
#9005-20

PAGE 65



**New!**


**PediFix Poster**  
Show Patients  
the Products  
You Recommend  
#LWALL1



**New!**

**Lamb's Wool by the Pound**  
Reliable Padding Patients Love  
#7970

PAGE 77



**New!**

**Visco-GEL® Rolled Digital Corn Pad**  
Targeted Gel Dot, Soft Ribbed Fabric  
#1010

PAGE 25



**New!**

**ClotrimaZoil®**  
New Package Look, Same Exclusive,  
Proven Formula!  
#3490

PAGE 59



**New!**

**DiabeticDefense™ Orthotic Insoles**  
Soft Support & Skin Protection.  
#2340

PAGE 48



**New!**

**Action Orthotics™**  
3/4-Length Affordable Support  
#P5610

PAGE 50

**Exclusive New Products**



**New!**  
**PressureOFF® Customizable Offloading Insoles**  
 Targeted Relief of Plantar Pathologies in Everyday Footwear #2400

PAGE 2



**New!**  
**Springthotics™ Carbon Graphite Rocker Sole Toe Spring Plates**  
 Reduce Joint Motion to Relieve Pain #4864

PAGE 10



**New!**  
**1stRaythotics™ Preforms®**  
 Integrated Morton's Extension, Arch Support Limit Joint Motion to Relieve Pain #4870

PAGE 7



**New!**  
**Halluxthotics™**  
 Rigid Morton's Extension Arch Support #4868

PAGE 9



**New!**  
**1stRaythotics™ Hallux Trainer**  
 Morton's Extension Inside Full Foot Orthotic #4876

PAGE 8



**New!**  
**1/2" Lateral Sole Wedge Insoles**  
 Control Supination, Relieve Peroneal Tendon, Knee Pain #2301

PAGE 11



**New!**  
**Just For Toenails®**  
 Tea Tree Oil Enriched  
 O<sub>2</sub> Permeable, 10-Chemical Free, Healthy Nail Formula #M905

PAGE 60



**New!**  
**Customizables® Bunion Fix™**  
 Personalized Relief, Comfort & Shoe Fit #1315-CTZ

PAGE 28



**New!**  
**Customizables® ToeBuddy®**  
 Adjustable Between Toe Separation #1035-CTZ

PAGE 29



**New!**  
**Customizables® Toe Spreader**  
 Adjustable Toe Alignment & Separation #1030-CTZ

PAGE 28



**New!**  
**Customizables® Hammer ToeCrutch®**  
 Personalized Toe Tip Lift #1037-CTZ

PAGE 29



**New!**  
**AdjustaLoop™ Visco-GEL® Toe Crest**  
 Lifts Toe Tips, Prevents Band Overtightening #3037-CTZ

PAGE 29

## Since 1885

### 4 family generations...

Expect outstanding products, service and support from the only fourth-generation, family-owned business in the footcare industry.



Heinrich Berkemann  
Company Founder, 1885



Dennis Case, Jon Case, Christopher Case

## RECOVERY

That's our theme this year. The world, our country, our profession, and our patients — we're all recovering from the unprecedented conditions of the Covid pandemic. I know you join me in an optimistic view — we can and will conquer this challenge. Together. With the help of science and medicine.

RECOVERY is also an appropriate theme for 30 new products you'll find inside, including some of our most exciting and potentially beneficial innovations ever. And for a company that's helped many millions of patients for well over a century, I don't say that lightly.

We expect PressureOFF™ Customizable Offloading Insoles will save limbs and lives, pain and suffering for thousands more this year. So much so that we're introducing them a bit prematurely — we're still perfecting our designs — and I'm eager for your feedback to make them better.

They'll benefit your diabetic patients striving to prevent ulcerations, protect wounds in remission, and relieve common plantar pathologies you treat frequently. By offloading on both plantar and dorsal insole surfaces, pressure and friction are reduced dramatically in your patient's everyday footwear, so they'll be more compliant, with better outcomes.

Also an innovation in recovery, new 1stRaythotics™ address an underserved patient need — convenient, comfortable relief of a variety of MTPJ conditions.

The hallux is implicated in more ailments than any other structure in the human foot. But until now, treatment options were limited. These new prefabricated orthotics feature integrated midfoot support and Morton's extensions that stabilize first ray motion. See variations — the Metatarsal & Neuroma Offloader, Reverse Morton's Extension, and Neutral styles. Aerospace-grade Carbon Fiber construction means thin, light, firm performance and comfort.

See our complimentary assortment of Halluxthotics®, Springthotics® spring plates, and other motion stabilizing insoles that will keep your patients active, with less pain — from professional athletes to geriatrics and those in between. All help MTPJ ailments, and can play an especially important role post operatively, to support today's most popular surgical constructs.

Patient visits will rebound, as we adapt to evolving Covid concerns and strategies. Your practice will return to full capacity. And we'll be here to help. With great products, 4 family generations of our family's commitment to you, your practice, your patients and your profession.

Wishing you well on your recovery journey,

Jon P. Case  
President

The PediFix Footcare Company

P.S. Recommend PediFix with confidence. Dispense from your office, Prescribe through locally-owned retailers. Enable patient direct orders as needed. Just be sure your patients get 'the good stuff'.

## To order or contact us:

**Call:** 1-800-424-5561 or 845-278-3800

**Fax:** 1-800-431-7801 or 845-277-2851

**Email:** [info@pedifix.com](mailto:info@pedifix.com)

**Mail:** 301 Fields Lane,  
Brewster, NY 10509

**Contact:** Your favorite authorized PediFix® Distributor

**Visit:** <https://medical.pedifix.com>

To access exclusive products and wholesale prices, log in or create a Healthcare Professional account.

## Unconditional Timeless Guarantee

If for any reason, you or your patients are ever dissatisfied with the performance, function or comfort of any PediFix product, return it at any time for a replacement, credit or full refund.

PRSRT STD  
 U.S. POSTAGE  
 PAID  
 RIPON, WI  
 PERMIT NO. 100

Change Service Requested

PediFix Footcare Company  
 301 Fields Lane  
 Brewster, NY 10509 USA



**SoreSpot® Silicone Tape**  
*Foot Care, Skin Care, Scar Care, Wound Care, Anywhere!*  
 #3095

PAGE 12

**AdjustaLoop™ Visco-GEL® Hammer Toe Crest**  
*Lifts sore toe tips, helps prevent band overtightening*  
 #3037

PAGE 29

**New!**  
**PressureOFF® Customizable Offloading Insoles**  
*Relieve Plantar Pressure, Friction & Pain with the First 'Removable Pegs' Offloading Insole System for Everyday Footwear*  
 #2402

PAGE 4



Call: 1-800-424-5561 • Visit: [medical.pedifix.com](http://medical.pedifix.com)  
 Email: [info@pedifix.com](mailto:info@pedifix.com) • Fax: 1-800-431-7801  
 Distributors Nationwide